

DIGISTIM III, Part# 8-1053-00

SPECIFICATIONS	
Double Burst	
Twitch	1 pulse/sec, 2 pulse/sec or 1 pulse/10 sec
Tetanus	50 Hz or 100 Hz
Train-of-Four	Single or automatic repeat every 10 sec
Output Current	adjustable 0-70 mA Hi Range, 2K Ohm Load, Adjustable 0-6 mA Lo Range
Simulus Pulse	
Power Source	One 9-volt alkaline battery
Size	2.09"H x 3.60"W x 5.75"D
Weight	Eleven (11) ounces including battery

DIGISTIM 3 PLUS, Part# 8-1053-02

SPECIFICATIONS	
Double Burst	2 60-ms bursts of 50 Hz separated by .75 s; repeats every 10 s
Twitch	1 s (1 pulse/s); 2 Hz (2 pulses/s)
Tetanus	50 Hz (50 pulses/s); 100 Hz (100 pulses/s)
Train-of-Four	4 pulses/2 s; repeats every 10 s
Output Current	0-70 mA, adjustable in 0.1-mA increments (HI range) 0-6.0 mA, adjustable in 0.1-mA increments (LO range)
Simulus Pulse	Square Eave Monophasic pulses
Power Source	One 9-volt alkaline battery
Size	2.09"H x 3.60"W x 5.75"D
Weight	Eleven (11) ounces including battery

DIGISTIM 2 PLUS, Part# 8-1053-72

SPECIFICATIONS	
Double Burst	2 60-ms bursts of 50 Hz separated by .75 s; single or automatic repeat every 10 s
Twitch	1 s (pulse/s); single or repeat
Tetanus	50 Hz (50 pulses/s); 100 Hz (100 pulses/s)
Train-of-Four	Single (4 pulses/2 s); repeat (4 pulses/2 s; repeats every 10 s)
Output Current	0-70 mA, adjustable in 0.1-mA increments (HI range) 0-6.0 mA, adjustable in 0.1-mA increments (LO range)
Simulus Pulse	Square Eave Monophasic pulses
Power Source	One 9-volt alkaline battery
Size	2.5"H x 6.06"W x 6.25"D
Weight	1 lb. 10 oz., including battery

MICROSTIM, part# 8-1053-60

SPECIFICATIONS	
Double Burst	
Twitch	1 pulse per second
Tetanus	100 Hz
Train-of-Four	4 Pulses per 2 Seconds
Output Current	Variable, 0-70 mA
Simulus Pulse	
Power Source	One 9-volt alkaline battery
Size	1.10"H x 2.42"W x 3.88"D
Weight	5.5 oz., including battery

MICROSTIM PLUS, part# 8-1053-62

SPECIFICATIONS	
Double Burst	2 60-ms bursts of 50 Hz seaprated by .75 s
Twitch	1 s (1 pulse/s)
Tetanus	50 Hz (50 pulses/s) or 100 Hz (100 pulses/s)
Train-of-Four	Single (4 pulses/2 s)
Output Current	Adjustable from 0-70 mA
Simulus Pulse	Square Eave Monophasic pulses
Power Source	One 9-volt alkaline battery
Size	1.10"H x 2.42"W x 3.88"D
Weight	5.5 oz., including battery