

Laryngeal
Mask Airways

Pages 34-37



Filtration


Pages 46-47

ET Tubes

Pages 30-33



INSIDE: Shop For All Your Anesthesia Care Products



End Tidal Cannulas

Pages 42-45


Nerve
Stimulators

Pages 18-19



Nasopharyngeal
Airways

Pages 24-25



**ON
SALE!**
HEPA & Filters
See Page 3



FEATURED PRODUCT

Capnography Solutions



ETCO₂ /O₂ CANNULA + REFLECTIVE CONNECTOR

- Dual nasal cannula simultaneously provides oxygen to both nostrils while obtaining CO₂ sampling during spontaneous breathing
- Designed to prevent mixing of fresh oxygen with CO₂
- Hydroscopic filter helps minimize moisture from being pulled into the monitor at the sampling port
- 3-Channel tubing prevents kinking

ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
5107F-SE	ETCO ₂ / O ₂ Nasal Cannula	Adult	7'	25	\$172.95
5110F-SE	ETCO ₂ / O ₂ Nasal Cannula	Adult	10'	25	\$190.95
5114F-SE	ETCO ₂ / O ₂ Nasal Cannula	Adult	14'	25	\$208.95



DIVIDED END TIDAL CANNULA

- True ETCO₂ sampling and provided O₂ delivery without limitations at higher O₂ rates
- Soft tubing for patient comfort
- Divided design separates gases
- Salter Eyes maintains sampling even during a cannula occlusion

ITEM	DESCRIPTION	LENGTH	PK	\$
4950FT-7-7-25	Cannula with Male Luer Connector	7'	25	\$49.95
4950FT-TG7-7-25	Cannula with Threaded Grip Connector	7'	25	\$50.95



CAPNOGRAPHY CANNULA

- Designed with minimal dead space in sampling line; resulting in a clear, accurate capnogram
- Comfort Soft Plus[®] tubing delivers premium patient comfort, reducing pressure ulcers and promoting patient wellness

ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
0912	CO ₂ /O ₂ Nasal Cannula	Adult	10'	25	\$246.95
0914	CO ₂ /O ₂ Nasal Cannula	Adult	14'	25	\$258.95
0916	CO ₂ /O ₂ Nasal Cannula	Pediatric	10'	25	\$235.95

INSIDE

Laryngoscopes 4-16
 Diagnostics 17-23
 Airway Management 24-41
 Respiratory 42-57
 Intubation 58-62
 Surgical 63-65
 Central Sterile 66-67



To place an order:
 Just give us a call - **1-800-237-5481**
 Send by fax - **1-800-671-7988**
 Via email - **info@bay-medical.com**
 Order online at **www.bay-medical.com**



Payment options:
 Credit Cards accepted: Visa, MasterCard, American Express, and Discover.

If you would like to set up an account, call us. Doctors, Nurses, Therapists and other Healthcare Professionals receive the same discounted prices given to Hospitals and Surgery Centers.

Terms and Conditions:
 If you receive items that do not meet your expectations, you may return them if in salable condition to Bay Medical for credit or replacement. Please contact our Customer Care Department at 1-800-237-5481 for a return authorization.

Full credit will be given for items returned in salable condition within 30 days. Special orders and sterile items are not returnable. Restocking fees may apply. Freight charges are not refundable.

We make every effort maintain the prices printed in this catalog and to provide you with the most up-to-date and accurate information. However, occasionally we receive price increases from our suppliers. Therefore prices are subject to change without notice. Please refer to our website for the most current pricing.

Contact information:
 Bay Medical
 2710 Northridge Dr. NW, Suite B
 Grand Rapids, MI 49544 USA

Toll Free: **1-800-237-5481**
 Fax: **1-800-671-7988**
 Email: **info@bay-medical.com**

Filtration Solutions for Oral Surgery



- ## B/V FILTER
- 99.99% Filtration efficacy
 - 34 mL Deadspace
 - > 250 mL VT range
 - 22 mm ID x 15 mm ID/22 mm OD



ITEM	DESCRIPTION	PK	REG	NOW ONLY \$24 YOU SAVE \$26
FH603003	Bacterial / Viral Filter	50	\$49.95	



- ## B/V FILTER & ADAPTER
- 99.99% Filtration efficacy
 - 34 mL Deadspace
 - > 250 mL VT range
 - 22 mm ID x 15 mm ID/22 mm OD



ITEM	DESCRIPTION	PK	REG	NOW ONLY \$38 YOU SAVE \$32
FH603003-CFF	Bacterial / Viral Filter & Adapter	50	\$69.95	



- ## B/V FILTER & TUBING
- 99.99% Filtration efficacy
 - 34 mL Deadspace
 - > 250 mL VT range
 - 22 mm ID x 15 mm ID/22 mm OD



ITEM	DESCRIPTION	PK	REG	NOW ONLY \$49 YOU SAVE \$51
FH603003-CFT	Bacterial / Viral Filter & Tubing	50	\$99.95	



- ## HEPA FILTER
- Effectively filters for COVID-19
 - 99.9999% Filtration efficacy
 - Sampling port
 - 62 mL Deadspace
 - 250 mL – 1500 mL Tidal volume
 - 22 mm OD/15 mm ID x 22 mm ID/15 mm OD



ITEM	DESCRIPTION	PK	REG	NOW ONLY \$95 YOU SAVE \$75
FH603074	HEPA / COVID-19 Filter	50	\$169.95	



- ## HEPA FILTER
- 22 mm ID x 15 mm ID / 22 mm OD
 - 99.99% Filtration efficacy
 - 40.6 mL Deadspace
 - >150 mL Tidal volume



ITEM	DESCRIPTION	PK	REG	NOW ONLY \$80 YOU SAVE \$90
FH603026	HEPA Filter	50	\$169.95	



- ## HEPA FILTER
- 99.9999% Filtration efficacy
 - Compatible with all brands of circuits
 - CO₂ port
 - Universal connectors
 - 33.58 mL Deadspace
 - 150 mL – 1500 mL Vt range



ITEM	DESCRIPTION	PK	REG	NOW ONLY \$70 YOU SAVE \$88
6159	HEPA Filter & Port	50	\$157.95	

FEATURED PRODUCT

SunOne™ SYSTEM

One time. One choice. One solution.

SINGLE PATIENT USE, FULLY DISPOSABLE, STERILE BLADES AND HANDLES.



SunMed

GREENLINE® / D™ LED HANDLE

- One complete, fully-contained unit
- Textured surface ensures secure grip
- Mercury-free, pre-installed batteries



ITEM	DESCRIPTION	PK	\$
5-5338-09	Medium Textured Grip Disposable Handle	20	\$153.95
5-5338-10	Penlite Textured Grip Disposable Handle	20	\$153.95
5-5338-11	Stubby Textured Grip Disposable Handle	20	\$153.95
5-5338-12	Standard Textured Grip Disposable Handle	20	\$153.95
5-0236-69	Large Ribbed Grip Disposable Handle (Clinically Clean)	20	\$197.95



ON CONTRACT

Premier PP-OR-1424
Vizient MS5172
CSF-LeeSar • R01-RMS02155
CALL FOR A QUOTE ON THIS AND OTHER PRODUCTS ON CONTRACT

SunMed

GREENLINE® / D™ MACINTOSH BLADE

- Non-painted, lightweight stainless steel
- Prevents cross contamination
- Polished acrylic stem produces exceptional illumination
- Green System compatible
- Perfect for hospital or EMS



ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
5-5332-01	Macintosh, Infant	1	92mm	20	\$104.95
5-5332-02	Macintosh, Child	2	100mm	20	\$104.95
5-5332-03	Macintosh, Medium Adult	3	130mm	20	\$104.95
5-5332-35	Macintosh, Extended Medium Adult	3.5	144mm	20	\$104.95
5-5332-04	Macintosh, Large Adult	4	155mm	20	\$104.95



\$5.24ea

SunMed

GREENLINE® / D™ MILLER BLADE

- Non-painted, lightweight stainless steel
- Prevents cross contamination
- Polished acrylic stem produces exceptional illumination
- Green System compatible
- Perfect for hospital or EMS



ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
5-5333-20	Miller, Premature	00	67mm	20	\$104.95
5-5333-00	Miller, Neonate	0	77mm	20	\$104.95
5-5333-01	Miller, Infant	1	102mm	20	\$104.95
5-5333-02	Miller, Child	2	153mm	20	\$104.95
5-5333-03	Miller, Medium Adult	3	195mm	20	\$104.95
5-5333-04	Miller, Large Adult	4	205mm	20	\$104.95



\$5.24ea

SunMed

GREENLINE® / D™ FLEX TIP BLADE

- Precise control of lever elevates epiglottis
- Non-pinching tip design
- Prevents cross contamination
- Green system compatible
- Great in urgent care, ambulances



ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
5-5382-03	Flex Tip, Medium Adult	3	195mm	10	\$197.95
5-5382-04	Flex Tip, Large Adult	4	205mm	10	\$197.95



\$19.79ea




SunMed
DISPOSABLE FIBER OPTIC LARYNGOSCOPE KITS


- Variety of blade and handle configurations
- All blades are single patient use and fully disposable
- Components in original, sterile packaging inside sturdy canvas pouch




ITEM	DESCRIPTION	PK	\$
5-5332-55	2 Disp Std Handles; Mac Blades 2, 3, 4; Miller Blades 2, 3	1	\$60.95
5-5332-57	1 Disp Std Handle; Disp Stubby Handle; Mac Blades 2, 3, 4; Miller Blades 2, 3	1	\$60.95
5-5336-54	2 Disp Penlite Handles; All-Metal Mac Blade 1; All Metal Miller Blades 00, 0, 1	1	\$60.95
5-5333-58	2 Disp Std Handles; Mac Blades 2, 3, 4; Miller Blades 0, 1, 2, 3	1	\$73.95
5-5333-57	1 Chrome Plated Med Handle; Mac Blades 2, 3, 4; Miller Blades 0, 1, 2, 3	1	\$104.95
5-5335-47	1 Chrome Plated Med Handle; All Metal Mac Blades 2, 3, 4; All Metal Miller Blades 0, 1, 2, 3	1	\$142.95
5-5133-57	LED Conventional Kit: 1 Med Handle; Mac English Blades: 2, 3, 4; Miller English Blades: 0, 1, 2, 3	1	\$98.95




5-5332-55




5-5332-57



5-5336-54



5-5333-58



5-5333-57



5-5133-57 (Conventional Laryngoscopes)

SunMed
GREENLINE® / D™ ALL-METAL DISPOSABLE MACINTOSH/MILLER BLADES

- Help reduce cross-contamination
- Eliminates reprocessing costs of blades
- SunOne System Component
- Fiber optic bundle protected in black plastic sheath
- Designed with three ball bearings in the heel for secure handle attachment
- Macintosh blades feature beaded tip reduces tissue trauma
- Constructed of surgical grade 303/304 non-painted stainless steel



ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
5-5335-00	All-Metal Macintosh, Neonate	0	80mm	20	\$129.95
5-5335-01	All-Metal Macintosh, Infant	1	88.5mm	20	\$129.95
5-5335-02	All-Metal Macintosh, Child	2	109.5mm	20	\$129.95
5-5335-03	All-Metal Macintosh, Medium Adult	3	131mm	20	\$129.95
5-5335-35	All-Metal Macintosh, Extended Med. Adult	3.5	145mm	20	\$129.95
5-5335-04	All-Metal Macintosh, Large Adult	4	160mm	20	\$129.95



5-5339-02

\$6.49ea



5-5339-01

SunMed
GREENLINE® / D™ ALL-METAL DISPOSABLE PHILLIPS BLADES

- Helps reduce cross-contamination
- Eliminates reprocessing costs of blades
- SunOne System Component
- Fiber optic bundle protected in black plastic sheath
- Designed with three ball bearings in the heel for secure handle attachment
- Constructed of surgical grade 303/304 non-painted stainless steel



ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
5-5339-01	All-Metal Phillips, Infant	1	113mm	20	\$129.95
5-5339-02	All-Metal Phillips, Child	2	166mm	20	\$129.95




\$6.49ea

ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
5-5336-20	All-Metal Miller, Premature	00	65mm	20	\$129.95
5-5336-00	All-Metal Miller, Neonate	0	78.5mm	20	\$129.95
5-5336-01	All-Metal Miller, Infant	1	104.5mm	20	\$129.95
5-5336-15	All-Metal Miller, Small Child	1.5	130mm	20	\$129.95
5-5336-02	All-Metal Miller, Child	2	157.3mm	20	\$129.95
5-5336-03	All-Metal Miller, Med. Adult	3	196.5mm	20	\$129.95
5-5336-04	All-Metal Miller, Large Adult	4	208mm	20	\$129.95



SunMed
GREENLINE® / D™ ALL-METAL DISPOSABLE WIS-HIPPLE BLADES

- Helps reduce cross-contamination
- Eliminates reprocessing costs of blades
- SunOne System Component
- Flexible fiber optic bundle protected in black plastic sheath
- Designed with three ball bearings in the heel for secure handle attachment
- Constructed of surgical grade 303/304 non-painted stainless steel



ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
5-5334-15	All-Metal Wis-Hipple, Small Child	1.5	115mm	20	\$129.95

\$6.49ea

LARYNGOSCOPES



PREMIER PRO DISPOSABLE BLADES

- Green System compatible
- Constructed of surgical stainless steel to reduce cleaning and cross-contamination
- ISO 7376 Compliant



MILLER BLADES

ITEM	DESCRIPTION	SIZE	PK	\$
62500	Premier Pro Miller	00	20	\$66.95
6250	Premier Pro Miller	0	20	\$66.95
6251	Premier Pro Miller	1	20	\$66.95
6252	Premier Pro Miller	2	20	\$66.95
6253	Premier Pro Miller	3	20	\$66.95
6254	Premier Pro Miller	4	20	\$66.95

MACINTOSH BLADES

ITEM	DESCRIPTION	SIZE	PK	\$
6241	Premier Pro Macintosh	1	20	\$66.95
6242	Premier Pro Macintosh	2	20	\$66.95
6243	Premier Pro Macintosh	3	20	\$66.95
62435	Premier Pro Macintosh	3.5	20	\$66.95
6244	Premier Pro Macintosh	4	20	\$66.95



PREMIER PRO DISPOSABLE HANDLES

- Green System compatible
- LED illumination
- Textured surface for better grip
- ISO 7376 Compliant



ITEM	DESCRIPTION	SIZE	PK	\$
6290	Premier Pro LED Handle	Stubby	20	\$113.95
6291	Premier Pro LED Handle	Penlite	20	\$113.95
6292	Premier Pro LED Handle	Medium	20	\$113.95



SunMed TOOTH PROTECTOR

- Latex free thermoplastic construction
- Helps avoid trauma to teeth



ITEM	DESCRIPTION	PK	\$
9-3041-03	Adult Tooth Protector	50	\$79.95



DUAL LED LARYNGOSCOPES

- Ultraviolet combined with white LED light provides clear illumination and definition exceptional depth perception and detail
- Illumination in the blade, powered by the handle
- Outstanding visualization of glottic opening
- Single patient use blades and handles help reduce risk of cross-contamination
- Eliminate reprocessing costs and save significantly
- Medical grade stainless steel blades and aircraft aluminum handles deliver quality that clinicians demand
- Full range of styles and sizes
- Medically clean



\$6.70ea

LARYNGOSCOPES

MACINTOSH DISPOSABLE BLADES

ITEM	DESCRIPTION	SIZE	PK	\$
1020.C	Stainless Steel Macintosh, Neonate	0	20	\$119.95
1021.C	Stainless Steel Macintosh, Infant	1	20	\$119.95
1022.C	Stainless Steel Macintosh, Child	2	20	\$119.95
1023.C	Stainless Steel Macintosh, Medium Adult	3	20	\$119.95
1029.C	Stainless Steel Macintosh, Extended Medium Adult	3	20	\$119.95
1024.C	Stainless Steel Macintosh, Large Adult	4	20	\$119.95

PHILLIPS DISPOSABLE BLADES

ITEM	DESCRIPTION	SIZE	PK	\$
2001.C	Stainless Steel Phillips, Infant	1	20	\$119.95
2002.C	Stainless Steel Phillips, Child	2	20	\$119.95


ROBERTSHAW DISPOSABLE BLADES

ITEM	DESCRIPTION	SIZE	PK	\$
2010.C	Stainless Steel Robertshaw, Neonate	0	20	\$119.95
2011.C	Stainless Steel Robertshaw, Infant	1	20	\$119.95

DISPOSABLE HANDLES

ITEM	DESCRIPTION	SIZE	PK	\$
2015.C	LED Aluminum Handle	Medium	20	\$133.95
2016.C	LED Aluminum Handle	Pediatric	20	\$133.95
2017.C	LED Aluminum Handle	Stubby	20	\$133.95
1058	Storage Bag	-	1	\$14.95

Significantly Improved Visualization of Airway with IntuBrite:



Competitive Fiber Optic LED

IntuBrite Combined Ultraviolet + White LED

\$5.99ea



MILLER DISPOSABLE BLADES

ITEM	DESCRIPTION	SIZE	PK	\$
1040.C	Stainless Steel Miller, Premature	000	20	\$119.95
1039.C	Stainless Steel Miller, Small Neonate	00	20	\$119.95
1030.C	Stainless Steel Miller, Neonate	0	20	\$119.95
1031.C	Stainless Steel Miller, Infant	1	20	\$119.95
1035.C	Stainless Steel Miller, Small Child	1.5	20	\$119.95
1032.C	Stainless Steel Miller, Child	2	20	\$119.95
1033.C	Stainless Steel Miller, Medium Adult	3	20	\$119.95
1034.C	Stainless Steel Miller, Large Adult	4	20	\$119.95

LARYNGOSCOPES



ON CONTRACT
Premier PP-OR-1424
Vizient MSS172
CSF-LeeSar
ROI-RMS02155
CALL FOR A QUOTE ON THIS AND OTHER
PRODUCTS ON CONTRACT



GREENLINE® MACINTOSH AMERICAN PROFILE BLADES

- Lightweight stainless steel
- No electrical contacts in the blade
- Fiber optic bundle provides brighter, cool white light
- Green System compatible
- Economical set includes popular sizes, medium handle, lamp, case



ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
5-5231-00	Macintosh American, Neonate	0	80mm	1	\$50.95
5-5231-01	Macintosh American, Infant	1	87mm	1	\$50.95
5-5231-02	Macintosh American, Child	2	108mm	1	\$50.95
5-5231-03	Macintosh American, Medium Adult	3	128mm	1	\$50.95
5-5231-35	Macintosh American, Adult	3.5	144mm	1	\$50.95
5-5231-04	Macintosh American, Large Adult	4	159mm	1	\$50.95
5-5231-44	Macintosh American Set: Medium Handle, Blade Sizes 1, 2, 3, 4, Lamp, Case				\$296.95



GREENLINE® MACINTOSH ENGLISH PROFILE BLADES

- Lightweight stainless steel
- No electrical contacts in the blade
- Fiber optic bundle provides brighter, cool white light
- Green System compatible
- Economical set includes popular sizes, medium handle, lamp and case



ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
5-5232-01	Macintosh English, Infant	1	92mm	1	\$65.95
5-5232-02	Macintosh English, Child	2	100mm	1	\$65.95
5-5232-03	Macintosh English, Medium Adult	3	130mm	1	\$65.95
5-5232-04	Macintosh English, Adult	4	155mm	1	\$65.95
5-5232-05	Macintosh English, Large Adult	5	175mm	1	\$65.95
5-5243-03	Macintosh English, Left Handed	3	130mm	1	\$65.95
5-5232-44	Macintosh English Set: Medium Handle, Blade Sizes 1, 2, 3, 4, Lamp, Case				\$296.95



GREENLINE® MACINTOSH GERMAN PROFILE BLADES

- Low profile
- Lightweight stainless steel
- No electrical contacts in the blade
- Green System compatible
- No cracks or crevices helps prevent cross contamination
- Large 5.0 mm integral fiber optic bundle provides more light
- Economical set includes popular sizes, medium handle, lamp, case

ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
5-5272-01	Macintosh German, Infant	1	93mm	1	\$82.95
5-5272-02	Macintosh German, Child	2	114mm	1	\$82.95
5-5272-03	Macintosh German, Medium Adult	3	135mm	1	\$82.95
5-5272-04	Macintosh German, Large Adult	4	155mm	1	\$82.95
5-5272-13	Strong Curve Macintosh German Profile	3	135mm	1	\$82.95
5-5272-44	Macintosh German Profile Set: Medium Handle, Blade Sizes 1, 2, 3, 4, Lamp, Case				\$395.95



EMERALD™ MACINTOSH GERMAN PROFILE BLADES

- 5.0mm bundle provides brightest light
- Matte finish reduces glare
- Integrated bundles design makes for easy cleaning; no crevices or hard to reach areas
- Green System compatible
- Lifetime warranty

ITEM	DESCRIPTION	SIZE	PK	\$
5-5211-01	Emerald Macintosh German, Infant	1	1	\$117.95
5-5211-02	Emerald Macintosh German, Child	2	1	\$117.95
5-5211-03	Emerald Macintosh German, Medium Adult	3	1	\$117.95
5-5211-35	Emerald Macintosh German, Extended Large Adult	3.5	1	\$117.95
5-5211-04	Emerald Macintosh German, Large Adult	4	1	\$117.95

LARYNGOSCOPES



MACINTOSH GERMAN PROFILE BLADES

- One-piece stainless steel design eliminates bacterial entrapment
- Sealed fiber optic bundle
- Manufactured to be sterilized Steam, STERRAD®, STERIS® or EO
- Green System compatible
- Pediatric blades feature wider blade design

ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
5-1991-01	Macintosh German, Infant	1	93mm	1	\$179.95
5-1991-02	Macintosh German, Child	2	115mm	1	\$174.95
5-1991-03	Macintosh German, Medium Adult	3	135mm	1	\$174.95
5-1991-04	Macintosh German, Large Adult	4	160mm	1	\$180.95
5-1991-05	Macintosh German, Extra Large Adult	5	176mm	1	\$215.95



SunMed
GREENLINE®
MILLER ENGLISH PROFILE BLADES

- Lightweight stainless steel
- Size of flange is reduced to minimize trauma
- Fiber optic light bundle provides brighter, cool white light
- Green System compatible
- Economical set includes popular sizes, medium handle, lamp, case

ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
5-5233-20	Miller English, Premature	00	67mm	1	\$65.95
5-5233-00	Miller English, Neonate	0	77mm	1	\$65.95
5-5233-01	Miller English, Infant	1	102mm	1	\$65.95
5-5233-15	Miller English, Small Child	1.5	130mm	1	\$65.95
5-5233-02	Miller English, Child	2	153mm	1	\$65.95
5-5233-03	Miller English, Medium Adult	3	195mm	1	\$65.95
5-5233-04	Miller English, Large Adult	4	205mm	1	\$65.95
5-5233-45	Miller English Profile Set: Medium Handle, Blades Sizes 0, 1, 2, 3, 4, Lamp, Case				\$329.95



SunMed
GREENLINE®
MILLER AMERICAN PROFILE BLADES

- Light on right, as viewed by user
- Lightweight stainless steel
- Size of flange is reduced to minimize trauma
- Fiber optic light bundle provides brighter, cool white light
- Green System compatible
- Replacement fiber optic bundle available

ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
5-5230-02	Miller American, Child	2	155mm	1	\$50.95
5-5230-03	Miller American, Medium Adult	3	195mm	1	\$50.95



propper
MILLER ENGLISH PROFILE BLADES

- Integrated fiber optic streamlined design
- No external light tubes for ease of cleaning
- Designed to be sterilized using steam, STERRAD®, STERIS® or EO
- Green System compatible
- Pediatric blades feature wider blade design

ITEM	DESCRIPTION	SIZE	PK	\$
5-1992-00	Miller English, Neonate	0	1	\$179.95
5-1992-01	Miller English, Infant	1	1	\$179.95
5-1992-15	Miller English, Small Child	1.5	1	\$179.95
5-1992-02	Miller English, Child	2	1	\$179.95
5-1992-03	Miller English, Medium Adult	3	1	\$179.95
5-1992-04	Miller English, Large Adult	4	1	\$179.95

SunMed
GREENLINE® MILLER GERMAN PROFILE BLADES

- Large integral fiber optic bundle provides more light
- No electrical contacts in the blade
- No cracks or crevices helps prevent cross contamination
- Green System compatible
- Economical set includes popular sizes, medium handle, lamp, case

ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
5-5273-20	Miller German, Premature	00	65mm	1	\$82.95
5-5273-00	Miller German, Neonate	0	77mm	1	\$82.95
5-5273-01	Miller German, Infant	1	110mm	1	\$82.95
5-5273-15	Miller German, Pediatric	1.5	130mm	1	\$82.95
5-5273-02	Miller German, Child	2	152mm	1	\$82.95
5-5273-03	Miller German, Medium Adult	3	192mm	1	\$82.95
5-5273-04	Miller German, Large Adult	4	207mm	1	\$82.95
5-5273-45	Miller German Profile Set: Medium Handle, Blade Sizes 0, 1, 2, 3, 4, Lamp, Case				\$461.95



RUSCH
EMERALD™ MILLER BLADES

- 5.0mm bundle provides brightest light
- Matte finish reduces glare
- Green System compatible
- Lifetime warranty
- Integrated bundles design makes for easy cleaning – no crevices or hard to reach areas

ITEM	DESCRIPTION	SIZE	PK	\$
5-5210-20	Emerald Miller, Small Premature	00	1	\$117.95
5-5210-00	Emerald Miller, Premature	0	1	\$117.95
5-5210-01	Emerald Miller, Infant	1	1	\$117.95
5-5210-02	Emerald Miller, Child	2	1	\$117.95
5-5210-03	Emerald Miller, Medium Adult	3	1	\$117.95
5-5210-04	Emerald Miller, Large Adult	4	1	\$117.95



RUSCH
GREEN SPEC® HANDLES

- Strong, durable stainless steel construction
- Bright Xenon lamp
- Textured surface to enhance grip
- Compatible with Rusch Emerald blades



ITEM	DESCRIPTION	SIZE	PK	\$
5-0242-61	Green Spec Handle	Penlite	1	\$199.95
5-0242-62	Green Spec Handle	Stubby	1	\$199.95
5-0242-63	Green Spec Handle	Medium	1	\$199.95

LARYNGOSCOPES

LARYNGOSCOPES



SUPERLUME LED HANDLES

- High-intensity LED handles provide a whiter, brighter light
- Lasts 1,000 times longer than halogen
- Superlume 12,000 hour LED light source eliminates need for replacement bulbs
- Available in sizes C, AA, and Stubby
- Batteries included

ITEM	DESCRIPTION	BATTERIES	PK	\$
5-1995-83	Penlite LED Fiber Optic Handle	2-AA	1	\$289.95
5-1995-82	Medium LED Fiber Optic Handle	2-C	1	\$273.95
5-1995-81	Stubby LED Fiber Optic Handle	2-AA	1	\$251.95



5-1995-83 5-1995-82 5-1995-81



INTEGRATED FIBER OPTIC HANDLES

- Unique removable battery sleeve for ease of cleaning and sterilization
- Compatible with all Green System fiber optic blades
- 5 Year warranty
- Batteries included



5-1995-01 5-1995-03 5-1995-02

ITEM	DESCRIPTION	BATTERIES	PK	\$
5-1995-01	Penlite Fiber Optic Handle	2-AA	1	\$256.95
5-1995-03	Medium Fiber Optic Handle	2-C	1	\$192.95
5-1995-02	Stubby Fiber Optic Handle	2-AA	1	\$264.95
5-0238-60	2.5V Xenon Reflector Bulb	-	1	\$16.95



5-0236-99 5-0236-90 5-0236-91

SunMed SUNBRITE GREENLINE® LED HANDLES

- 35,000 hours LED life
- Chrome-plated brass
- Batteries not included



5-0240-62 5-0240-63

ITEM	DESCRIPTION	BATTERIES	PK	\$
5-0236-99	Medium LED Fiber Optic Handle	2-C	1	\$109.95
5-0236-90	Penlite LED Fiber Optic Handle	2-AA	1	\$109.95
5-0236-91	Stubby LED Fiber Optic Handle	2-AAA	1	\$109.95
5-0240-62	F/O ECO LED Bulb Cartridge, Orange	-	1	\$23.95
5-0240-63	F/O Halogen/Xenon Bulb Cartridge, Yellow	-	1	\$51.95



GREENLINE® CHROME-PLATED/STAINLESS STEEL HANDLES

- Compatible with all Green System blades
- Xenon halogen reflector lamp produces the whitest brightest illumination
- Batteries not included



5-0240-61

CHROME PLATED BRASS

ITEM	DESCRIPTION	BATTERIES	PK	\$
5-0236-11	Stubby Fiber Optic Handle	2-AA	1	\$64.95
5-0236-10	Penlite Fiber Optic Handle	2-AA	1	\$64.95
5-0236-09	Medium Fiber Optic Handle	2-C	1	\$64.95
5-0236-15	Large Fiber Optic Handle	2-D	1	\$64.95

STAINLESS STEEL WITH SATIN FINISH

ITEM	DESCRIPTION	BATTERIES	PK	\$
5-0236-36	Micro Lithium Fiber Optic Handle	Lithium 123	1	\$75.95
5-0236-38	Mini Lithium Fiber Optic Handle	Lithium 123	1	\$75.95
5-0236-30	Penlite Fiber Optic Handle	2-AA	1	\$71.95
5-0236-31	Stubby Fiber Optic Handle	2-AA	1	\$71.95
5-0236-39	Medium Fiber Optic Handle	2-C	1	\$75.95
5-0236-35	Large Fiber Optic Handle	2-D	1	\$75.95
5-0240-61	F/O Halogen/Xenon Bulb Cartridge, Blue	-	1	\$20.95



FIBER OPTIC HOWLAND LOCK

- Assists with making difficult intubations possible
- Built in leverage prevents prying and reduces the possibility of broken teeth
- Fits all Green System laryngoscope handles and blades

ITEM	DESCRIPTION	PK	\$
5-0236-16	Fiber Optic Howland Lock	1	\$60.95



GREENLINE® PREMIUM LED HANDLES

- Autoclavable without having to remove lamp
- Twice as bright as Sunbrite LED handles
- Waterproof
- Chrome-plated brass with satin finish
- Long battery life
- Lifetime warranty

ITEM	DESCRIPTION	BATTERIES	PK	\$
5-0236-80	Penlite Premium LED Handle	2-AA	1	\$142.95
5-0236-89	Medium Premium LED Handle	2-C	1	\$142.95
5-0240-63	F/O Halogen/Xenon Bulb Cartridge, Yellow	-	1	\$51.95





Stop packing your own airway kits.
We have them—ready to use!

SunMed
COMPLETE INTUBATION KITS

- Complete intubation kits available with conventional, fiber optic or disposable laryngoscope options
- All components in original sterile packaging

ITEM	DESCRIPTION	PK	\$
8-1050-85	Conventional Intubation Kit	1	\$373.95
8-1050-95	Fiber Optic Intubation Kit	1	\$736.95
8-1050-94	Disposable Fiber Optic Blade Intubation Kit	1	\$296.95
8-1050-97	Disposable Fiber Optic Handle & Blade Intubation Kit	1	\$219.95
8-1050-80	Nylon Case Only	1	\$98.95

KIT INCLUDES:

- 8 Color Coded Guedel Airways (sizes 40mm-100mm)
- 4 Nasopharyngeal Airways (sizes 20, 24, 28, 32FR)
- 2 Uncuffed ET Tubes (sizes 2.5, 4mm)
- 8 Cuffed ET Tubes (sizes 5-10mm)
- 1 ET Tube Introducer (Bougie) with Coudé tip
- 1 Adult Stylette
- 1 Adult Magill Forceps
- 1 Pediatric Magill Forceps
- 4 Packets of Sterile Lubricating Jelly
- 1 Syringe (10ml)
- 1 Nylon Case
- 2 "AA" & "C" Batteries (if applicable)
- Replacement laryngoscope lamps (if applicable)

LARYNGOSCOPE HANDLES & BLADES:

- 1 Penlite Laryngoscope Handle
- 1 Medium Laryngoscope Handle
- 5 Miller Laryngoscope Blades (sizes 0-4)
- 4 Macintosh Laryngoscope Blades (sizes 1-4)

LARYNGOSCOPE COMBINATIONS + KIT ITEM CODES



SunMed
CRNA CLINICAL TRAINING KIT

- Essential diagnostic tools for anesthesia professionals in training
- Includes nerve stimulator, monoscope and four chest pieces with double adhesive disks
- Excellent for residency work and beyond at a reasonable student price
- One convenient purchase rather than separate product orders
- Durable canvas pouch keeps all tools close at hand

ITEM	DESCRIPTION	PK	\$
6-1030-50	Basic Training Kit Includes SunStim Nerve Stimulator	1	\$245.95
6-1030-52	Student Training Kit Includes SunStim Plus Nerve Stimulator	1	\$256.95
6-1030-53	Professional Training Kit Includes SunStim Pro Nerve Stimulator	1	\$295.95

WITH THE CRNA KIT YOU SAVE \$22

EACH KIT CONTAINS:

1 Peripheral Nerve Stimulator
Your choice of:



1 Disposable Monoscope



4 Chest Pieces (with 8 double adhesive disks)



Professional CRNA Training Kit
6-1030-53



Basic CRNA Training Kit
6-1030-50



Student CRNA Training Kit
6-1030-52

DIAGNOSTICS

DIAGNOSTICS



Leading Brand of Nerve Stimulators

Brand New with Digital Screen!



SUNSTIM™ ORIGINAL, SINGLE TETANUS STIMULATOR

SUNSTIM™ PLUS DOUBLE BURST, DUAL TETANUS FUNCTIONALITY

SUNSTIM™ PRO LCD DIGITAL SCREEN ERGONOMIC, COMFORTABLE SHAPE COLORFUL, PROTECTIVE SILICONE WRAP

SUNSTIM™

- Constant voltage stimulator for monitoring neuromuscular block
- Evaluates muscle relaxant dosage during surgery and recovery

SUNSTIM™ PLUS

- Complete muscle relaxant that fits in the palm of your hand
- Delivers sufficient output to ensure supramaximal stimulation; output current is fully adjustable from 0-70mA
- Features 50 and 100Hz tetanus, automatic train-of-four, and double-burst stimulation

SUNSTIM™ PRO

- Ergonomic and compact shape; fits better in hand and more comfortable to use
- Easy-to-read digital screen display for DBS, Tetanus, TOF and Twitch functionality
- Removable silicone wrap for secure grip; helps protect device
- Compatible with all SunStim accessories
- Includes gray silicone wrap



8-1053-GN GREEN

8-1053-TL TEAL

SILICONE PRO WRAP™

- Protects device
- Provides a secure grip
- Color code by clinician or floor/wing
- Personalize - available in 9 colors

ITEM	DESCRIPTION	PK	\$
See item numbers at left	Silicone Pro Wrap™	1	\$21.95



Made exclusively for the SunStim™ Pro



SunStim 8-1053-60

SunStim Plus 8-1053-62

SunStim Pro 8-1053-63

SPECIFICATIONS

FEATURES	SUNSTIM™ - 8-1053-60	SUNSTIM™ PLUS - 8-1053-62	SUNSTIM™ PRO - 8-1053-63
DBS	N/A	Two (2) 60 ms bursts of 50 Hz separated by 0.75 secs (750 milliseconds)	Two (2) 60 ms bursts of 50 Hz separated by 0.75 secs (750 milliseconds)
Tetanus	100 Hz	50 Hz to 100 Hz	50 Hz to 100 Hz
TOF	4 pulses in a period of 2 seconds	4 pulses in a period of 2 seconds	4 pulses in a period of 2 seconds
Output Current	Adjustable from 0-70mA by control knob	Adjustable from 0-70mA by control knob	Adjustable from 0-70mA by +/- switch
Power Source	One 9V alkaline battery	One 9V alkaline battery	One 9V alkaline battery
Functionality Indicators	Pulse and Battery illumination	Pulse and Battery illumination	LCD display screen
Pro Wrap™	N/A	N/A	Silicone Pro Wrap - available in 9 colors



8-1053-87



8-1051-78

8-1051-79



8-1053-76

8-1053-85

STIMULATORS

ITEM	DESCRIPTION	PK	\$
8-1053-60	SunStim™ Peripheral Nerve Stimulator	1	\$235.95
8-1053-62	SunStim™ Plus Peripheral Nerve Stimulator	1	\$246.95
8-1053-63	SunStim™ Pro Peripheral Nerve Stimulator	1	\$285.95

ACCESSORIES

ITEM	DESCRIPTION	PK	\$
8-1053-87	Zippered Pouch	1	\$10.95
8-1053-76	Stim Ball Electrodes & Hollow Pin	1	\$32.95
8-1051-78	6' Stim Lead Snap & Shielded Pin	1 pr	\$38.95
8-1051-79	6' Stim Lead Clip & Shielded Pin	1 pr	\$38.95
8-1053-85	Pole Mount Bracket for 8-1053-60 & 8-1053-62	1	\$144.95
8-1053-95	Pole Mount Bracket for 8-1053-63 (not shown)	1	\$144.95
8-1053-81	Stim Overlay for 8-1053-60 (not shown)	1	\$10.95
8-1053-82	Stim Overlay for 8-1053-62 (not shown)	1	\$10.95

All nerve stimulators come ready to use with an instruction manual, battery, ball electrodes and 6' leads with clips.

They all include a 1 year limited warranty; call for details!

DIAGNOSTICS



Call for exact product configurations, monitor compatibility and great pricing!
1-800-237-5481



COMFORT CHECK BLOOD PRESSURE CUFFS

- Soft white fabric material for enhanced comfort
- Disposable – reduces costly disinfection procedures and contamination concerns
- Increase infection control with single patient use cuffs
- Meets AAMI and AHA Cuff Sizing Recommendations
- Color coded printing on each cuff for easy identification
- Full range of connectors and adapters allow standardization of cuffs
- Easy connection to electronic NiBP monitors or manual blood pressure measurement systems



SINGLE TUBE CONNECTOR OPTIONS	CONNECTOR	1-Tube Bayonet	1-Tube Screw Type	1-Tube Male Twist Lock	1-Tube Female Twist Lock	1-Tube Poppa
DOUBLE TUBE CONNECTOR OPTIONS	CONNECTOR	2-Tube Screw Type	2-Tube Male/Female Twist Lock	2-Tube Female Twist Lock	2-Tube Manual Inflation and Male Luer	2-Tube Poppa and Manual Inflation

CUFF SIZE	CUFF COLOR CODING	LIMB CIRCUMFERENCE
Infant	Brown	9-15 cm
Child	Green	14-21.5 cm
Small Adult	Light Blue	20.5-28.5 cm
Adult	Navy Blue	27-35 cm
Adult Long	Navy Blue	27-35 cm
Large Adult	Plum	34-44 cm
Large Adult Long	Plum	34-44 cm
Thigh	Brown	42-52 cm




ANEROID SPHYGMOMANOMETER KIT

- EZView™ luminescent dial
- Complete with standard air release valve and zippered case
- 5 year warranty-neoprene inflation system
- 10 year calibration warranty



ITEM	DESCRIPTION	PK	\$
2-9141-01	discontinued	1	\$30.95
2-9141-02	Child Sphygmomanometer	1	\$36.95
2-9141-03	Adult Sphygmomanometer	1	\$36.95
2-9141-04	Large Adult Sphygmomanometer	1	\$39.95
2-9141-05	Thigh Sphygmomanometer	1	\$43.95

DIAGNOSTICS



Only \$1.49ea




LC ACCU-BAR TEMPERATURE TREND INDICATOR

- Provides continuous, accurate temperature trending non-invasively
- Core adjusted 94°-105°F (35°- 40.5°C)
- Accuracy ± 1°F / 0.5°C



ITEM	DESCRIPTION	PK	\$
8-6505-20	LC Accu-Bar Temperature Trend Indicator	100	\$149.95



Only \$0.69ea




TEMP-ASSURE TEMPERATURE MONITOR

- The industry's only single-degree increment thermometers that feature the clearly defined normal temperature highlighted by the unique arrow
- Accurate and inexpensive
- Convenient dispenser box mounts anywhere
- 94°-104° F (35°-40°C)



ITEM	DESCRIPTION	PK	\$
8-6515-40	Temp-Assure Monitor	100	\$69.95



Only \$0.58ea




BROADBAND LIQUID CRYSTAL INDICATOR

- Widest range and yet economical; consolidate your use of multiple indicators with one solution
- Accurate and inexpensive
- 90°-106° F (33°-41°C)



ITEM	DESCRIPTION	PK	\$
8-6515-30	BroadBand Liquid Crystal Indicator	100	\$58.95



Only \$0.92ea



CRYSTALLINE® II TEMPERATURE INDICATOR

- Core adjusted 92°-106°F (33°- 41°C)
- Handy dispenser box mounts to wall or equipment



ITEM	DESCRIPTION	PK	\$
8-6501-01	Crystalline® II Temperature Indicator	100	\$92.95



Only \$2.06ea



STAT SKAN® LCD

- Core adjusted 95°-107°F (36°-41°C)
- Unique circular shape allows for precise localized monitoring



ITEM	DESCRIPTION	PK	\$
8-6510-10	Stat Skan® LCD Temp Indicator	100	\$206.95



Only \$0.48ea



MINI TEMP INDICATOR

- Clearly marked normal core temperatures
- Latex free, accurate, and inexpensive
- 96°-104° F (36°-40°C)



ITEM	DESCRIPTION	PK	\$
8-6515-20	Mini Temp Liquid Crystal Indicator	100	\$48.95



Only \$1.60ea



KWIK-SKAN™

- Core adjusted 92°-106°F (34°- 41°C)
- Oval design for more accurate readings
- Convenient dispensing box mounts easily to wall or anesthesia carts



ITEM	DESCRIPTION	PK	\$
8-6510-50	Kwik-Skan™ Temp Monitor	100	\$160.95

DIAGNOSTICS



TruMed Medical, Inc.

TYMPANIC PROBE

- Less-invasive alternative to esophageal or rectal probes for monitoring core body temperature
- Soft, memory-cell foam conforms to ear canal
- Seals out ambient air which can affect temperature reading



ITEM	DESCRIPTION	PK	\$
6-1038-64	Tympanic Probe 400 Series Sensor	100	\$206.95



TruMed Medical, Inc.

ESOPHAGEAL STETHOSCOPE

- Monitor temperature, heart and breathing of anesthetized patients during surgery
- Soft, thin-walled cuff provides maximum sound transmission
- Male luer allows convenient attachment to standard acoustical earpieces



ITEM	DESCRIPTION	OD	PK	\$
6-1037-09	Esophageal 400 Series Sensor / Standard Cuff	9Fr	50	\$209.95
6-1037-12	Esophageal 400 Series Sensor / Standard Cuff	12Fr	50	\$221.95
6-1037-18	Esophageal 400 Series Sensor / Standard Cuff	18Fr	50	\$196.95
6-1037-24	Esophageal 400 Series Sensor / Standard Cuff	24Fr	50	\$248.95




TruMed Medical, Inc.

SKIN SENSOR

- Accurately assess skin temperature – reflective metal cover insulates sensor to prevent OR lights or heat lamps from affecting actual temperature
- Foam pad seals area between skin and sensor to eliminate ambient air that can affect temperature reading
- Hypo-allergenic adhesive foam allows for easy placement and removal



ITEM	DESCRIPTION	PK	\$
6-1038-62	Skin Temp 400 Series Sensor	200	\$586.95



TruMed Medical, Inc.

400 SERIES MONITOR INTERFACE CABLES

- Compatible with all probes featured here – and many other standard probes – reducing the need for multiple cables
- Connects to most monitors
- 15' length for ease of use
- Reusable

ITEM	DESCRIPTION	PK	\$
6-1038-41	400 Series with HP Type Connector	1	\$35.95
6-1038-42	400 Series with Phonojack Connector	1	\$35.95

DIAGNOSTICS

Sensoronics

OXYGEN SENSOR REPLACEMENTS

- For most anesthesia and respiratory applications
- Unconditional 1 year warranty

ITEM	OXYGEN SENSOR	REPLACES	PK	\$
8-3102-51	SS-25A	Ceramatec CAG-25E/250E (Maxtec Max-25E/250E)	1	\$69.95
8-3101-41	SS-14A	Drager 68-50645	1	\$97.95
8-3101-42	SS-10A	Ohmeda 0237-2034-700	1	\$97.95
8-3102-01	SS-20A	Ohmeda 6050-0004-110	1	\$97.95
8-3101-11	SS-11A	Teledyne R17, Catalyst (MSA) 406931	1	\$69.95
8-3101-21	SS-12A	Teledyne R22, Catalyst (MSA) 472062	1	\$69.95
8-3100-10	SS-01	Teledyne T1 & T2	1	\$64.95
8-3100-30	SS-03	Teledyne T7	1	\$69.95

Detailed cross reference chart available. Call us!



SS-03
8-3100-30



SS-10A
8-3101-42



SS-11A
8-3101-11



SS-20A
8-3102-01



SS-14A
8-3101-41



SS-12A
8-3101-21



SS-01
8-3100-10

SunMed CHEST PIECES

- Components can be purchased separately
- Deep sound chambers
- Easily connect to ear pieces
- Lightweight stainless steel - newborn and infant sizes
- Stainless steel - child and adult sizes
- Chest pieces are 100% autoclavable
- Set includes: 1 adult chest piece, 1 child chest piece, 1 infant chest piece, 8 double adhesive disk and 1 Comfort FIT SoftEar™ monoscope



ITEM	DESCRIPTION	PK	\$
6-0025-00	7/8" Dia. Newborn Chest Piece	1	\$9.95
6-0025-01	1-1/8" Dia. Infant Chest Piece	1	\$9.95
6-0025-02	1-1/8" Dia. Child Chest Piece	1	\$10.95
6-0025-03	1-1/2" Dia. Adult Chest Piece	1	\$10.95
6-0025-25	Double Adhesive Disks	100/Roll	\$15.95
6-1030-25	Comfort Fit Soft Ear™ Monoscope	25	\$65.95
6-0025-05	Chest Piece Set	1	\$32.95





SunMed
LATEX ROBERTAZZI AIRWAY

- Smooth surface for easy Introduction
- Soft and flexible to reduce trauma
- Multi size kit available

SunMed
LATEX FREE AIRWAY

- Soft with adjustable flange for patient safety

SALTER LABS
NASO WITH SLIDING DEPTH RING

- Patented tip design facilitates insertion and tube advancement
- Flex-tips are clinically proven to reduce nasal bleeding
- Soft, flexible PVC with trumpet-style design

ITEM	SIZE	ID	OD	LENGTH	BOX OF 10
1-5076-20	20Fr	5.0mm	6.7mm	115mm	\$30.95
1-5076-22	22Fr	5.5mm	7.3mm	120mm	\$30.95
1-5076-24	24Fr	6.0mm	7.6mm	130mm	\$30.95
1-5076-26	26Fr	6.5mm	8.6mm	135mm	\$30.95
1-5076-28	28Fr	7.0mm	9.3mm	135mm	\$30.95
1-5076-30	30Fr	7.5mm	10.0mm	140mm	\$30.95
1-5076-32	32Fr	8.0mm	10.7mm	150mm	\$30.95
1-5076-34	34Fr	8.5mm	11.2mm	160mm	\$30.95
1-5076-36	36Fr	9.0mm	12.0mm	175mm	\$30.95
1-5076-04	4 pk-sizes 20, 24, 28 & 32 plus Sterile Lube. Jelly				\$14.95
1-5076-06	6 pk-sizes 22, 24, 26, 28, 30 & 32 plus Sterile Lube. Jelly				\$20.95
1-5076-09	9 pk-sizes 20, 22, 24, 26, 28, 30, 32, 34 & 36 plus Sterile Lube. Jelly				\$30.95

ITEM	SIZE	ID	OD	LENGTH	BOX OF 10
1-5072-12	12Fr	3.0mm	4.0mm	95mm	\$50.95
1-5072-14	14Fr	3.5mm	4.7mm	95mm	\$50.95
1-5072-16	16Fr	4.0mm	5.3mm	95mm	\$50.95
1-5072-18	18Fr	4.5mm	6.0mm	95mm	\$50.95
1-5072-20	20Fr	5.0mm	6.7mm	115mm	\$50.95
1-5072-22	22Fr	5.5mm	7.3mm	120mm	\$50.95
1-5072-24	24Fr	6.0mm	7.6mm	120mm	\$50.95
1-5072-26	26Fr	6.5mm	8.6mm	130mm	\$50.95
1-5072-28	28Fr	7.0mm	9.3mm	170mm	\$50.95
1-5072-30	30Fr	7.5mm	10.0mm	170mm	\$50.95
1-5072-32	32Fr	8.0mm	10.7mm	170mm	\$50.95
1-5072-34	34Fr	8.5mm	11.2mm	170mm	\$50.95
1-5072-36	36Fr	9.0mm	12.0mm	170mm	\$50.95

ITEM	OD	ID	BOX OF 10
H-NPA-12-10	12Fr	4.0mm	\$37.95
H-NPA-14-10	14Fr	4.7mm	\$37.95
H-NPA-16-10	16Fr	5.3mm	\$37.95
H-NPA-18-10	18Fr	6.0mm	\$37.95
H-NPA-20-10	20Fr	6.7mm	\$37.95
H-NPA-22-10	22Fr	7.3mm	\$37.95
H-NPA-24-10	24Fr	8.0mm	\$37.95
H-NPA-26-10	26Fr	8.7mm	\$37.95
H-NPA-28-10	28Fr	9.3mm	\$37.95
H-NPA-30-10	30Fr	10.0mm	\$37.95
H-NPA-32-10	32Fr	10.7mm	\$37.95
H-NPA-34-10	34Fr	11.3mm	\$37.95
H-NPA-36-10	36Fr	12.0mm	\$37.95

FREE Robertazzi Nasopharyngeal Airway provide a patent airway

zi Nasopharyngeal Airway provide a patent airway

al Airway provide a patent airway:

airway with water soluble lubricant. airway into nostril. Gently rotate the airway to ease entry.


28Fr
 ID: 7.0mm
 OD: 9.3mm
 Length: 135mm

30Fr
 ID: 7.5mm
 OD: 10.0mm
 Length: 135mm

32Fr
 ID: 8.0mm
 OD: 10.7mm
 Length: 145mm

Extra large size denoted on airway packaging for quick visual reference and accurate selection

Airway specifications easily identified to ensure proper NPA selection



SunMed
LATEX FREE ROBERTAZZI AIRWAY

- Flexible neoprene performs like latex airways
- Soft material prevents damage to nasal passage
- Blunt tip
- Individually wrapped

ITEM	SIZE	ID	OD	LENGTH	BOX OF 10
1-5075-20	20Fr	5.0mm	6.7mm	115mm	\$40.95
1-5075-22	22Fr	5.5mm	7.3mm	120mm	\$40.95
1-5075-24	24Fr	6.0mm	7.6mm	120mm	\$40.95
1-5075-26	26Fr	6.5mm	8.6mm	130mm	\$40.95
1-5075-28	28Fr	7.0mm	9.3mm	135mm	\$40.95
1-5075-30	30Fr	7.5mm	10.0mm	135mm	\$40.95
1-5075-32	32Fr	8.0mm	10.7mm	145mm	\$40.95
1-5075-34	34Fr	8.5mm	11.2mm	160mm	\$40.95
1-5075-36	36Fr	9.0mm	12.0mm	175mm	\$40.95

ITEM	DESCRIPTION	\$
1-5075-04	4 pk - sizes: 20, 24, 28 & 32 includes Sterile Lubricating Jelly Packs	\$21.95
1-5075-06	6 pk - sizes: 22, 24, 26, 28, 30 & 32 includes Sterile Lubricating Jelly Packs	\$30.95
1-5075-09	9 pk - sizes: 20, 22, 24, 26, 28, 30, 32, 34 & 36 includes Sterile Lubricating Jelly Packs	\$43.95

RUSCH
 A TELEFLEX COMPANY
ROBERTAZZI AIRWAY

- Soft-no damage to nasal passage
- Resists collapse or kinking
- Easy open blister package

ITEM	SIZE	ID	OD	LENGTH	BOX OF 10
1-5071-20	20Fr	5.0mm	6.7mm	100mm	\$48.95
1-5071-22	22Fr	5.5mm	7.3mm	105mm	\$48.95
1-5071-24	24Fr	6.0mm	7.6mm	110mm	\$48.95
1-5071-26	26Fr	6.5mm	8.6mm	118mm	\$48.95
1-5071-28	28Fr	7.0mm	9.3mm	125mm	\$48.95
1-5071-30	30Fr	7.5mm	10.0mm	127mm	\$48.95
1-5071-32	32Fr	8.0mm	10.7mm	120mm	\$48.95
1-5071-34	34Fr	8.5mm	11.2mm	155mm	\$48.95
1-5071-36	36Fr	9.0mm	12.0mm	170mm	\$48.95



AIRWAY MANAGEMENT



AIRWAY MANAGEMENT

SunMed
CLEAR PVC GUEDEL AIRWAY

- Soft, smooth flexible design for patient comfort
- Rigid bite block
- Individually wrapped
- Available as a set; ideal for crash carts and EMS



BOX OF 10	BOX OF 50	SIZE	COLOR	BOX OF 10	BOX OF 50
1-5010-30	1-1521-30	30mm / 3cm	Clear	\$14.95	\$49.95
1-5010-20	1-1521-20	40mm / 4cm	Pink	\$14.95	\$49.95
1-5010-10	1-1521-10	50mm / 5cm	Blue	\$14.95	\$49.95
1-5010-01	1-1521-01	60mm / 6cm	Black	\$14.95	\$49.95
1-5010-02	1-1521-02	70mm / 7cm	White	\$14.95	\$49.95
1-5010-03	1-1521-03	80mm / 8cm	Green	\$14.95	\$49.95
1-5010-04	1-1521-04	90mm / 9cm	Yellow	\$14.95	\$49.95
1-5010-05	1-1521-05	100mm / 10cm	Red	\$14.95	\$49.95
1-5010-06	1-1521-06	110mm / 11cm	Orange	\$14.95	\$49.95
1-5010-98	Complete PVC Guedel Airway Set, one of each size				\$14.95

PORTEX
SOFT CLEAR PVC AIRWAY

- Single use
- Smooth rounded edges for patient comfort
- Rigid nylon bite-block facilitates passage of a suction catheter



ITEM	mm	SIZE	COLOR	PK	\$
1-5012-40	2.67	4cm	Pink	10	\$15.95
1-5012-50	2.67	5cm	Dark Blue	10	\$15.95
1-5012-60	3.33	6cm	Black	10	\$15.95
1-5012-70	3.33	7cm	White	10	\$15.95
1-5012-80	4.0	8cm	Green	10	\$15.95
1-5012-90	4.0	9cm	Yellow	10	\$15.95
1-5012-10	4.0	10cm	Red	10	\$15.95
1-5012-11	4.67	11cm	Light Blue	10	\$15.95
1-5012-12	4.67	12cm	Orange	10	\$15.95

RUSCH
SOFT CLEAR PVC AIRWAY

- Flexible
- Individually wrapped



ITEM	SIZE	COLOR	PK	\$
1-5008-30	30mm / 3cm	Clear	10	\$20.95
1-5008-40	40mm / 4cm	Pink	10	\$20.95
1-5008-50	50mm / 5cm	Blue	10	\$20.95
1-5008-60	60mm / 6cm	Black	10	\$20.95
1-5008-70	70mm / 7cm	White	10	\$20.95
1-5008-80	80mm / 8cm	Green	10	\$20.95
1-5008-90	90mm / 9cm	Yellow	10	\$20.95
1-5008-99	100mm / 10cm	Red	10	\$20.95
1-5008-11	120mm / 12cm	Purple	10	\$20.95

SunMed
GUEDEL AIRWAY WITH BITE BLOCK

- Rigid bite block
- Individually wrapped
- Available as a set; ideal for crash carts and EMS



BOX OF 10	BOX OF 50	SIZE	COLOR	BOX OF 10	BOX OF 50
1-1504-40	1-1505-40	40mm / 4cm	Pink	\$7.95	\$30.95
1-1504-50	1-1505-50	50mm / 5cm	Blue	\$7.95	\$30.95
1-1504-60	1-1505-60	60mm / 6cm	Black	\$7.95	\$30.95
1-1504-70	1-1505-70	70mm / 7cm	White	\$7.95	\$30.95
1-1504-80	1-1505-80	80mm / 8cm	Green	\$7.95	\$30.95
1-1504-90	1-1505-90	90mm / 9cm	Yellow	\$7.95	\$30.95
1-1504-99	1-1505-99	100mm / 10cm	Red	\$7.95	\$30.95
1-1504-11	1-1505-11	110mm / 11cm	Orange	\$7.95	\$30.95
1-1504-00	Complete Guedel Airway with Bite Block Set, one of each size				\$7.95



CATH-GUIDE GUEDEL AIRWAY

- Flexible vinyl with rigid bite block
- Designed with three internal channels
- Individually wrapped



ITEM	SIZE	PK	\$
1-1160-55	55mm / 550cm	48	\$57.95
1-1160-60	60mm / 600cm	48	\$57.95
1-1160-70	70mm / 700cm	48	\$57.95
1-1160-80	80mm / 800cm	48	\$57.95
1-1160-90	90mm / 900cm	48	\$57.95
1-1160-10	100mm / 1000cm	48	\$57.95
1-1160-11	110mm / 1100cm	48	\$57.95
1-1160-12	120mm / 1200cm	48	\$57.95



BERMAN ORAL AIRWAY

- Full range of sizes available with length imprinted on airway
- Premium upgrade without a premium price
- Non-sterile



ITEM	SIZE	COLOR	PK	\$
1-5009-40	40mm / 4cm	Pink	10	\$5.95
1-5009-50	50mm / 5cm	Blue	10	\$5.95
1-5009-60	60mm / 6cm	Black	10	\$5.95
1-5009-70	70mm / 7cm	White	10	\$5.95
1-5009-80	80mm / 8cm	Green	10	\$5.95
1-5009-90	90mm / 9cm	Yellow	10	\$5.95
1-5009-10	100mm / 10cm	Red	10	\$5.95
1-5009-11	110mm / 11cm		10	\$5.95



FREE SHIPPING
for ONLINE ORDERS
of **\$500+**
www.bay-medical.com

Some restrictions apply.



VENTED SEPTUM BERMAN AIRWAY

- Designed to allow air flow from either side
- Available as a set; ideal for crash carts and EMS



BOX OF 10	BOX OF 50	Size	Color	BOX OF 10	BOX OF 50
1-1508-40	1-1506-40	40mm / 4cm	Pink	\$4.95	\$15.95
1-1508-50	1-1506-50	50mm / 5cm	Blue	\$4.95	\$15.95
1-1508-60	1-1506-60	60mm / 6cm	Black	\$4.95	\$15.95
1-1508-70	1-1506-70	70mm / 7cm	White	\$4.95	\$15.95
1-1508-80	1-1506-80	80mm / 8cm	Green	\$4.95	\$15.95
1-1508-90	1-1506-90	90mm / 9cm	Yellow	\$4.95	\$15.95
1-1508-99	1-1506-99	100mm / 10cm	Red	\$4.95	\$15.95
1-1508-11	1-1506-11	110mm / 11cm	Orange	\$4.95	\$15.95
1-1508-00	Complete Berman Airway Set, ten of each size				\$25.95



OVASSAPIAN FIBER OPTIC INTUBATING AIRWAY

- Designed primarily for use with a fiber optic endoscope
- Designed to allow for removal without disconnecting or displacing positioned endotracheal tube
- Distal half is open posteriorly to provide fiberscope maneuverability
- Individually packaged - fits tube sizes 8.0mm - 9.00mm



ITEM	DESCRIPTION	PK	\$
1-5001-03	Ovassapian Fiberoptic Intubating Airway - Adult	5	\$37.95



WILLIAMS AIRWAY

- Effectively intubate the trachea
- Helps guide fiber optic bronchoscope placement



ITEM	DESCRIPTION	SIZE	PK	\$
1-1503-09	Airway Intubator - Adult Female	90mm / 9cm	10	\$38.95
1-1503-10	Airway Intubator - Adult Male	100mm / 10cm	10	\$38.95



AIRWAY MANAGEMENT



AIRWAY MANAGEMENT

SunMed CUFFED WITH STYLETTE

- Barrel cuff profile minimizes capillary constriction within the trachea
- Aluminum stylette
- Radiopaque strip embedded for safety
- Smooth tube tip to minimize trauma during intubation
- High volume, low pressure cuff provides excellent tracheal wall seal
- Murphy eye



ITEM	Fr	ID	BOX OF 10
1-7343-30	12	3.0mm	\$30.95
1-7343-35	14	3.5mm	\$30.95
1-7343-40	16	4.0mm	\$30.95
1-7343-45	18	4.5mm	\$30.95
1-7343-50	20	5.0mm	\$30.95
1-7343-55	22	5.5mm	\$30.95
1-7343-60	24	6.0mm	\$30.95
1-7343-65	26	6.5mm	\$30.95
1-7343-70	28	7.0mm	\$30.95
1-7343-75	30	7.5mm	\$30.95
1-7343-80	32	8.0mm	\$30.95
1-7343-85	34	8.5mm	\$30.95
1-7343-90	36	9.0mm	\$30.95
1-7343-95	38	9.5mm	\$30.95
1-7343-10	40	10.0mm	\$30.95

SunMed UNCUFFED

- Oral/nasal use
- Radiopaque strip embedded for safety
- Smooth tube tip to minimize trauma during intubation
- Standard 15mm connector
- Murphy eye



ITEM	Fr	ID	BOX OF 10
1-7330-20	8	2.0mm	\$10.95
1-7330-25	10	2.5mm	\$10.95
1-7330-30	12	3.0mm	\$10.95
1-7330-35	14	3.5mm	\$10.95
1-7330-40	16	4.0mm	\$10.95
1-7330-45	18	4.5mm	\$10.95
1-7330-50	20	5.0mm	\$10.95
1-7330-55	22	5.5mm	\$10.95
1-7330-60	24	6.0mm	\$10.95
1-7330-65	26	6.5mm	\$10.95
1-7330-70	28	7.0mm	\$10.95

SunMed CUFFED

- Barrel cuff profile minimizes capillary constriction within the trachea
- Radiopaque strip embedded for safety
- Smooth tube tip to minimize trauma during intubation
- High volume, low pressure cuff provides excellent tracheal wall seal
- Murphy eye



ITEM	Fr	ID	BOX OF 10
1-7333-25	10	2.5mm	\$16.95
1-7333-30	12	3.0mm	\$16.95
1-7333-35	14	3.5mm	\$16.95
1-7333-40	16	4.0mm	\$16.95
1-7333-45	18	4.5mm	\$16.95
1-7333-50	20	5.0mm	\$16.95
1-7333-55	22	5.5mm	\$16.95
1-7333-60	24	6.0mm	\$16.95
1-7333-65	26	6.5mm	\$16.95
1-7333-70	28	7.0mm	\$16.95
1-7333-75	30	7.5mm	\$16.95
1-7333-80	32	8.0mm	\$16.95
1-7333-85	34	8.5mm	\$16.95
1-7333-90	36	9.0mm	\$16.95
1-7333-95	38	9.5mm	\$16.95
1-7333-10	40	10.0mm	\$16.95

SunMed ORAL CUFFED (AGT) (RAE) TYPE

- Directs tube downward to rest on patient's chin
- Oral use
- Murphy eye
- Barrel cuff profile minimizes capillary constriction within the trachea
- High volume, low pressure cuff provides excellent tracheal wall seal
- Smooth tube tip to minimize trauma during intubation
- Standard 15mm connector



ITEM	Fr	ID	BOX OF 10
1-7321-40	16	4.0mm	\$54.95
1-7321-45	18	4.5mm	\$54.95
1-7321-50	20	5.0mm	\$54.95
1-7321-55	22	5.5mm	\$54.95
1-7321-60	24	6.0mm	\$54.95
1-7321-65	26	6.5mm	\$54.95
1-7321-70	28	7.0mm	\$54.95
1-7321-75	30	7.5mm	\$54.95
1-7321-80	32	8.0mm	\$54.95
1-7321-85	34	8.5mm	\$54.95
1-7321-90	36	9.0mm	\$54.95

SunMed ORAL UNCUFFED (AGT) (RAE) TYPE

- Directs tube downward to rest on patient's chin
- Oral use
- Murphy eye
- Smooth tube tip to minimize trauma during intubation
- Standard 15mm connector



ITEM	Fr	ID	BOX OF 10
1-7320-30	12	3.0mm	\$43.95
1-7320-35	14	3.5mm	\$43.95
1-7320-40	16	4.0mm	\$43.95
1-7320-45	18	4.5mm	\$43.95
1-7320-50	20	5.0mm	\$43.95
1-7320-55	22	5.5mm	\$43.95

SunMed NASAL CUFFED (AGT) (RAE) TYPE

- Directs tube over patient's forehead, avoiding any possible facial trauma
- Nasal use
- Murphy eye
- Barrel cuff profile minimizes capillary constriction within the trachea
- High volume, low pressure cuff provides excellent tracheal wall seal
- Smooth tube tip to minimize trauma during intubation
- Standard 15mm connector



ITEM	Fr	ID	BOX OF 10
1-7323-40	16	4.0mm	\$54.95
1-7323-45	18	4.5mm	\$54.95
1-7323-50	20	5.0mm	\$54.95
1-7323-55	22	5.5mm	\$54.95
1-7323-60	24	6.0mm	\$54.95
1-7323-65	26	6.5mm	\$54.95
1-7323-70	28	7.0mm	\$54.95
1-7323-75	30	7.5mm	\$54.95
1-7323-80	32	8.0mm	\$54.95

SunMed NASAL UNCUFFED (AGT) (RAE) TYPE

- Directs tube over patient's forehead, avoiding any possible facial trauma
- Nasal use
- Murphy eye
- Smooth tube tip to minimize trauma during intubation
- Standard 15mm connector




ITEM	Fr	ID	BOX OF 10
1-7322-30	12	3.0mm	\$40.95
1-7322-35	14	3.5mm	\$40.95
1-7322-40	16	4.0mm	\$40.95
1-7322-45	18	4.5mm	\$40.95
1-7322-50	20	5.0mm	\$40.95
1-7322-55	22	5.5mm	\$40.95

Exceptionally Easy Intubation

Flexible, rounded tip slides gently and quickly along airway anatomy

Tip hugs intubation guides, introducers and stylets

Two Murphy eyes adjacent to tip




FLEX-TIP ENDOTRACHEAL TUBES

- Designed to facilitate smooth, rapid, non-traumatic intubation
- Unique, patented design has a flexible, curved, centered, tapered distal tip
- Tapered tip is designed to flex and slide gently along the airway through the larynx, down the trachea on its curved backside without catching on airway tissue



See and Experience the Difference!



FLEX-TIP TUBE

Slides gently along airway

No gap between tip of Parker tube and fiber bundle

STANDARD ET TUBE

Gap between leading edge of standard ET tube and fiber bundle. This gap can catch on airway anatomy as the tube is 'railroaded' downward.



CUFFED

ITEM	ID	PK	\$
H-PFHV-40-10	4.0mm	10	\$20.95
H-PFHV-45-10	4.5mm	10	\$20.95
H-PFHV-50-10	5.0mm	10	\$20.95
H-PFHV-55-10	5.5mm	10	\$20.95
H-PFHV-60-10	6.0mm	10	\$20.95
H-PFHV-65-10	6.5mm	10	\$20.95
H-PFHV-70-10	7.0mm	10	\$20.95
H-PFHV-75-10	7.5mm	10	\$20.95
H-PFHV-80-10	8.0mm	10	\$20.95
H-PFHV-85-10	8.5mm	10	\$20.95
H-PFHV-90-10	9.0mm	10	\$20.95
H-PFHV-95-10	9.5mm	10	\$20.95

CUFFED WITH STYLET

ITEM	ID	PK	\$
H-PFST-40-10	4.0mm	10	\$38.95
H-PFST-45-10	4.5mm	10	\$38.95
H-PFST-50-10	5.0mm	10	\$38.95
H-PFST-55-10	5.5mm	10	\$38.95
H-PFST-60-10	6.0mm	10	\$38.95
H-PFST-65-10	6.5mm	10	\$38.95
H-PFST-70-10	7.0mm	10	\$38.95
H-PFST-75-10	7.5mm	10	\$38.95
H-PFST-80-10	8.0mm	10	\$38.95
H-PFST-85-10	8.5mm	10	\$38.95
H-PFST-90-10	9.0mm	10	\$38.95
H-PFST-95-10	9.5mm	10	\$38.95

UNCUFFED

ITEM	ID	PK	\$
H-PFUO-25-10	2.5mm	10	\$23.95
H-PFUO-30-10	3.0mm	10	\$23.95
H-PFUO-35-10	3.5mm	10	\$23.95
H-PFUO-40-10	4.0mm	10	\$23.95
H-PFUO-45-10	4.5mm	10	\$23.95
H-PFUO-50-10	5.0mm	10	\$23.95
H-PFUO-55-10	5.5mm	10	\$23.95
H-PFUO-60-10	6.0mm	10	\$23.95
H-PFUO-65-10	6.5mm	10	\$23.95
H-PFUO-70-10	7.0mm	10	\$23.95

ORAL CUFFED

ITEM	ID	BOX OF 10
H-PFOC-40-10	4.0mm	\$48.95
H-PFOC-45-10	4.5mm	\$48.95
H-PFOC-50-10	5.0mm	\$48.95
H-PFOC-55-10	5.5mm	\$48.95
H-PFOC-60-10	6.0mm	\$48.95
H-PFOC-65-10	6.5mm	\$48.95
H-PFOC-70-10	7.0mm	\$48.95
H-PFOC-75-10	7.5mm	\$48.95
H-PFOC-80-10	8.0mm	\$48.95
H-PFOC-85-10	8.5mm	\$48.95
H-PFOC-90-10	9.0mm	\$48.95
H-PFOC-95-10	9.5mm	\$48.95

ORAL UNCUFFED

ITEM	ID	BOX OF 10
H-PFOU-30-10	3.0mm	\$45.95
H-PFOU-35-10	3.5mm	\$45.95
H-PFOU-40-10	4.0mm	\$45.95
H-PFOU-45-10	4.5mm	\$45.95
H-PFOU-50-10	5.0mm	\$45.95
H-PFOU-55-10	5.5mm	\$45.95
H-PFOU-60-10	6.0mm	\$45.95
H-PFOU-65-10	6.5mm	\$45.95
H-PFOU-70-10	7.0mm	\$45.95
H-PFOU-75-10	7.5mm	\$45.95
H-PFOU-80-10	8.0mm	\$45.95
H-PFOU-85-10	8.5mm	\$45.95
H-PFOU-90-10	9.0mm	\$45.95

NASAL UNCUFFED

ITEM	ID	BOX OF 5
H-PFNU-30-5	3.0mm	\$34.95
H-PFNU-35-5	3.5mm	\$34.95
H-PFNU-40-5	4.0mm	\$34.95
H-PFNU-45-5	4.5mm	\$34.95
H-PFNU-50-5	5.0mm	\$34.95
H-PFNU-55-5	5.5mm	\$34.95
H-PFNU-60-5	6.0mm	\$34.95
H-PFNU-65-5	6.5mm	\$34.95
H-PFNU-70-5	7.0mm	\$34.95

NASAL CUFFED

ITEM	ID	BOX OF 5
H-PFNC-40-5	4.0mm	\$34.95
H-PFNC-45-5	4.5mm	\$34.95
H-PFNC-50-5	5.0mm	\$34.95
H-PFNC-55-5	5.5mm	\$34.95
H-PFNC-60-5	6.0mm	\$34.95
H-PFNC-65-5	6.5mm	\$34.95
H-PFNC-70-5	7.0mm	\$34.95
H-PFNC-75-5	7.5mm	\$34.95
H-PFNC-80-5	8.0mm	\$34.95
H-PFNC-85-5	8.5mm	\$34.95
H-PFNC-90-5	9.0mm	\$34.95



Helps Deliver Safety & Confidence in Airway Management



AIR-Q®3 INTUBATING LARYNGEAL AIRWAY

- Unique ET tube ramp and epiglottis elevator aid in fast, easy, and safe intubations
- Designed to optimize airway placement and seal
- Anatomically correct placement of the masked airway
- Lifts the epiglottis to maximize access
- Guide ramp assists with ET tube insertion and placement

GREAT PRICING AVAILABLE!
CALL 800-237-5481



AIR-Q®3 STANDARD, MANUAL CUFF INFLATION

ITEM	DESCRIPTION	SIZE	MAX OETT	TUBE LENGTH	INFLATION VOLUME	PK	\$
10-30105	Air-Q3 Infant, 4-7 kg	1.0	4.5 mm	9.0 cm	0.5 - 1.0 mL	10	CALL FOR PRICING
10-30155	Air-Q3 Infant/Pediatric, 7-17 kg	1.5	5.0 mm	11.0 cm	1.0 mL	10	CALL FOR PRICING
10-30205	Air-Q3 Pediatric, 17-30 kg	2.0	5.5 mm	14.0 cm	1.0 - 2.0 mL	10	CALL FOR PRICING
10-30305	Air-Q3 Adult, 30-60 kg	3.0	7.0 mm	16.0 cm	2.0 - 3.0 mL	10	CALL FOR PRICING
10-30405	Air-Q3 Adult, 60-80 kg	4.0	8.0 mm	18.0 cm	3.0 - 4.0 mL	10	CALL FOR PRICING
10-30505	Air-Q3 Adult, >80 kg	5.0	9.0 mm	20.0 cm	4.0 - 5.0 mL	10	CALL FOR PRICING

AIR-Q®3G STANDARD, MANUAL CUFF INFLATION, GASTRIC ACCESS


ITEM	DESCRIPTION	SIZE	MAX OETT	TUBE LENGTH	INFLATION VOLUME	PK	\$
10-50105	Air-Q3 Infant, 4-7 kg	1.0	4.5 mm	9.0 cm	0.5 - 1.0 mL	10	CALL FOR PRICING
10-50155	Air-Q3 Infant/Pediatric, 7-17 kg	1.5	5.0 mm	11.0 cm	1.0 mL	10	CALL FOR PRICING
10-50205	Air-Q3 Pediatric, 17-30 kg	2.0	5.5 mm	14.0 cm	1.0 - 2.0 mL	10	CALL FOR PRICING
10-50305	Air-Q3 Adult, 30-60 kg	3.0	7.0 mm	16.0 cm	2.0 - 3.0 mL	10	CALL FOR PRICING
10-50405	Air-Q3 Adult, 60-80 kg	4.0	8.0 mm	18.0 cm	3.0 - 4.0 mL	10	CALL FOR PRICING
10-50505	Air-Q3 Adult, >80 kg	5.0	9.0 mm	20.0 cm	4.0 - 5.0 mL	10	CALL FOR PRICING



AIR-Q®3



AIR-Q®3 PROVIDES BETTER VISUAL CLARITY



Competitive ILA

Air-Q3 Series with Epiglottis Elevator



AIR-Q®3

Provides exceptional ease of ETT placement when unexpected intubation is required.



AIR-Q®3G

Features a self-pressurizing cuff with an extra wide gastric inlet that accommodate NG tubes up to 18 Fr and two gastric channels.



AIR-Q®SP3*


Uses positive pressure ventilation to self-pressurize the cuff. Cuff pressures cycle between peak airway pressure and PEEP levels.



AIR-Q®SP3G*

Uses positive pressure ventilation to self-pressurize the cuff. Cuff pressures cycle between peak airway pressure and PEEP levels. Features a self-pressurizing cuff with an extra wide gastric inlet that accommodate NG tubes up to 18 Fr and two gastric channels.

Perfect Placement with Better Visual Clarity



Gastric access channels on both sides of airway tube in Air-Q3G series

Color-coded connector for easy size identification
Removable tethered connector enables insertion of endotracheal tube through airway tube.

Anatomically curved, medical grade silicone airway tube with integrated bite block

Aperture enables transfer of air if the inlet becomes blocked

Exclusive epiglottis elevator lifts the epiglottis to maximize cuff placement

Soft medical grade silicone cuff


Self-pressurizing cuff inflation is controlled by inspiratory/expiratory pressure. No manual inflation needed (for Air-QSP3* and Air-QSP3G* models only)

Extra large gastric inlet available in Air-Q3G series

Built-up mask heel for improved seal

ET tube ramp helps facilitate intubation and directs the endotracheal tube toward the laryngeal inlet

*This Air-Q3 series is available online: www.bay-medical.com




DISPOSABLE, COST-EFFECTIVE MASK AIRWAY

- Constructed of medical-grade PVC
- Engineered using extra soft matte textured materials to replicate the performance of silicone
- Reduces the risk of cross-contamination
- Easy to insert while providing a secure seal
- Pilot balloon indicates cuff pressure
- Color-coded pilot balloon makes for quick and accurate size identification



ITEM	DESCRIPTION	SIZE	WEIGHT	PK	\$
1-2062-10	LarySeal™ Clear PVC	1	< 11 lbs	10	\$115.95
1-2062-15	LarySeal™ Clear PVC	1.5	11-22 lbs	10	\$115.95
1-2062-20	LarySeal™ Clear PVC	2	22-44 lbs	10	\$115.95
1-2062-25	LarySeal™ Clear PVC	2.5	44-66 lbs	10	\$115.95
1-2062-30	LarySeal™ Clear PVC	3	66-110 lbs	10	\$115.95
1-2062-40	LarySeal™ Clear PVC	4	110-154 lbs	10	\$115.95
1-2062-50	LarySeal™ Clear PVC	5	154-220 lbs	10	\$115.95



SINGLE USE, SILICONE MASK AIRWAY

- Performs comparable to that of reusable silicone mask airways
- Single-use
- Reduces the risk of cross-contamination
- Easy to insert while providing a secure seal
- Pilot balloon indicates cuff pressure
- Color-coded pilot balloon makes for quick and accurate size identification



ITEM	DESCRIPTION	SIZE	WEIGHT	PK	\$
1-2061-10	LarySeal™ Blue Silicone	1	<11 lbs	10	\$171.95
1-2061-15	LarySeal™ Blue Silicone	1.5	11-22 lbs	10	\$171.95
1-2061-20	LarySeal™ Blue Silicone	2	22-44 lbs	10	\$171.95
1-2061-25	LarySeal™ Blue Silicone	2.5	44-66 lbs	10	\$171.95
1-2061-30	LarySeal™ Blue Silicone	3	66-110 lbs	10	\$171.95
1-2061-40	LarySeal™ Blue Silicone	4	110-154 lbs	10	\$171.95
1-2061-50	LarySeal™ Blue Silicone	5	154-220 lbs	10	\$171.95



AURASTRAIGHT™ AIRWAY

- Clearly marked guidelines for easy insertion, cuff inflation volume and patient weight on pilot balloon
- Maximum intracuff pressure 60cm H₂O for all sizes
- MRI compatible



Item	Description	size	Weight	MAX CUFF VOLUME	BOX OF 10
1-2252-01	Neonate/Infant	1	up to 5kg	4ml	\$106.95
1-2252-15	Infant	1.5	5kg - 10kg	7ml	\$106.95
1-2252-02	Pediatric	2	10kg - 20kg	10ml	\$106.95
1-2252-25	Pediatric	2.5	20kg - 30kg	14ml	\$106.95
1-2252-03	Adult	3	30kg - 50kg	20ml	\$106.95
1-2252-04	Adult	4	50kg - 70kg	30ml	\$106.95
1-2252-05	Adult	5	70kg - 100kg	40ml	\$106.95
1-2252-06	Adult	6	100kg and over	50ml	\$106.95



AURA40™ REUSABLE STRAIGHT AIRWAY

- Steam autoclavable version of Ambu's popular AuraOnce
- Reusable silicone-ease of use and cleaning
- Clearly marked guidelines for easy insertion
- Maximum intracuff pressure 60cm H₂O for all sizes



ITEM	DESCRIPTION	SIZE	MAX CUFF	PK	\$
1-2265-01	Neonate up to 5kg	1	4ml	1	\$240.95
1-2265-15	Infants 5kg - 10kg	1.5	7ml	1	\$240.95
1-2265-02	Infant/Child 10kg - 20kg	2	10ml	1	\$240.95
1-2265-25	Children 20kg - 30kg	2.5	14ml	1	\$240.95
1-2265-03	Children 30kg - 50kg	3	20ml	1	\$240.95
1-2265-04	Adults 50kg - 70kg	4	30ml	1	\$240.95
1-2265-05	Adults 70kg - 100kg	5	40ml	1	\$240.95

AIRWAY MANAGEMENT



RÜSCH
A TELEFLEX COMPANY

ENDOTEST CUFF PRESSURE MANOMETER

- Fills, monitors and adjusts cuff pressure of low pressure endotracheal, endobronchial and tracheostomy tubes
- Helps minimize damage to the mucous membranes of the tracheal wall when used to monitor cuff pressure
- Easy-to-read-gauge features a green zone that highlights the recommended inflation pressure range
- One-hand operation
- Fingertip trigger switch increases and maintains desired cuff pressure
- Designed to hang from nearby anesthesia equipment



VBM Medical Inc.

POCKET CUFF PRESSURE GAUGE

- For inflation and pressure control on tracheal tubes with high volume, low pressure cuffs
- Complete with 100 cm connecting tubing
- Adult range 0-120cm H₂O, green range 20-30cm H₂O
- Pediatric gauge has pressure range from 0-60 cm H₂O, recommended green range from 5 to 15cm H₂O and release valve for sensitive adjustment of the desired cuff pressure

8-0050-64

8-0050-60

ITEM	DESCRIPTION	PK	\$
8-0050-60	Pocket Cuff Pressure Gauge with Hook	1	\$197.95
8-0050-64	Universal Pressure Gauge	1	\$197.95

ITEM	DESCRIPTION	PK	\$
8-0050-00	EndoTest Cuff Pressure Manometer	1	\$320.95

Immediate Visualization of CO₂ Expiration



Adult R500P30

Suitable for use with masks or ET tubes during resuscitation



SunMed

MAXCAP CO₂ INDICATOR

- ETCO₂ and flow indication device designed for quick visualization of exhaled CO₂ and expiratory flow
- 360° visualization of color change which can be seen from any location around patient
- High definition color changing foam material (purple to yellow)
- Foam will not dry out unlike litmus paper used in many indicators
- Weighs less than other CO₂ indicators for less pull on ET tubes
- Available in 3 sizes



Infant R500P32



Pediatric R500P31

ITEM	DESCRIPTION	SIZE	WT	PK	\$
R500P30	MaxCap CO ₂ Indicator	Adult	>15 kg (33 lbs)	10	\$83.95
R500P31	MaxCap CO ₂ Indicator	Pediatric	1 kg - 15 kg	10	\$100.95
R500P32	MaxCap CO ₂ Indicator	Infant	0.25 kg - 6 kg	10	\$170.95


SunMed

IGUARD™ EYE PROTECTOR


- Clear plastic protective covers with non-allergenic self-adhesive foam
- Protects patients eyes during general anesthesia, against eye abrasion, lacerations and scleral hemorrhage
- Fast and accurate application
- Individually packaged



ITEM	DESCRIPTION	PK	\$
9-0210-00	Adult IGuard	25	\$146.95
9-0210-21	Child IGuard	25	\$148.95
9-0210-20	Adult Double Foam IGuard	25	\$148.95



Thicker adhesive foam is more comfortable for patient, even in prone position.



Tab design means one-handed release of adhesive paper.



AIRWAY MANAGEMENT



FACE IT™ FACE SHIELDS


- Protect your face!
- Meet's OSHA's blood born Pathogen Regulations
- Comfortable self shaping foam headpiece provides eye protection
- Crystal clear plastic with anti-fog coating ensures clear viewing
- Drape shield has fluid barrier fabric with hook and loop closure for secure coverage



ITEM	DESCRIPTION	PK	\$
9-0300-01	Full Face Shield	50	\$185.95
9-0300-02	Drape Shield	50	\$227.95



AIRWAY MANAGEMENT



SunMed
FLAT BREATHING BAGS

- Reusable
- Fabricated from the finest neoprene
- 22mm bushing



ITEM	DESCRIPTION	SIZE	PK	\$
4-1041-05	Flat Breathing Bag	1/2 Liter	1	\$14.95
4-1041-10	Flat Breathing Bag	1 Liter	1	\$15.95
4-1041-15	Flat Breathing Bag	1-1/2 Liter	1	\$15.95
4-1041-20	Flat Breathing Bag	2 Liter	1	\$16.95
4-1041-30	Flat Breathing Bag	3 Liter	1	\$16.95
4-1041-40	Flat Breathing Bag	4 Liter	1	\$17.95
4-1041-50	Flat Breathing Bag	5 Liter	1	\$17.95



HY-TAPE®

- The Original Pink Tape® that is a unique water-resistant and washable tape especially suited for use with ET tubes
- Resists moisture, conforms to body contours and adheres firmly to oily and hairy skin
- Will not break down over time, so skin remains intact promoting faster healing and reducing infection
- 36 Individually wrapped rolls, each 5 yds long
- Available in three widths: 1/2", 3/4" and 1"



ITEM	DESCRIPTION	PK	\$
9-0075-05	HY-Tape - 1/2" x 5 yds.	36 rolls	\$85.95
9-0075-75	HY-Tape - 3/4" x 5 Yds.	36 rolls	\$160.95
9-0075-10	HY-Tape - 1" x 5 Yds.	36 rolls	\$160.95



SunMed
CHRISTMAS TREE CONNECTOR

- Color coded for appropriate gas connection identification
- One-piece design has no loose parts
- Two-piece design allows rotation while tubing is attached

ITEM	DESCRIPTION	COLOR	PK	\$
8-2311-20	One-Piece Connector	Green	50	\$21.95
8-2311-50	Two-Piece Connector	Green	50	\$27.95
8-2311-26	One-Piece Connector	Yellow	50	\$21.95
8-2311-22	One-Piece Connector	White	1000	\$353.95

AIRWAY MANAGEMENT



STRAIGHT

SunMed
REUSABLE STAINLESS STEEL ENDOTRACHEAL TUBE CONNECTORS

- High quality stainless steel



Straight	Acute 60° No Port	90° No Port	Size	Pk	Straight	Acute 60° No Port	90° No Port
3-3100-40	3-3101-40	3-3120-40	4mm x 15mm	1	\$16.95	\$27.95	\$27.95
3-3100-50	3-3101-50	3-3120-50	5mm x 15mm	1	\$16.95	\$27.95	\$27.95
3-3100-60	3-3101-60	3-3120-60	6mm x 15mm	1	\$16.95	\$27.95	\$27.95
3-3100-70	3-3101-70	3-3120-70	7mm x 15mm	1	\$16.95	\$27.95	\$27.95
3-3100-80	3-3101-80	3-3120-80	8mm x 15mm	1	\$16.95	\$27.95	\$27.95
3-3100-90	3-3101-90	3-3120-90	9mm x 15mm	1	\$16.95	\$27.95	\$27.95



ACUTE 60° NO PORT



90° NO PORT

NOW ONLY \$15
3-3120-50,
3-3120-70

FREE SHIPPING
for ONLINE ORDERS
of \$500+
www.bay-medical.com

Some restrictions apply.

Ethox
BITE BLOCK

- Assists with naso-gastric tube placement
- Accommodates up to 40Fr tube
- Non-sterile
- Individually wrapped



ITEM	DESCRIPTION	PK	\$
2060	Bite Block	20	\$48.95



HUDSON RCI
BITEGUARD™ BITE BLOCK

- Rigid, yet compressible design prevents damage to LMAs, dental work and fingertips
- Convenient handle provides quick simple insertion



ITEM	DESCRIPTION	PK	\$
9-3042-40	Biteguard™ Bite Block	50	\$42.95

RESPIRATORY



As Low As
\$6.71ea



STAT ORANGE CAPNOGRAPHY CANNULA

- True ETCO₂ sampling with proven O₂ delivery
- Features reflective connector comparable to standard orange and long-term yellow end tidal cannulas
- Reflective connector compatible with leading capnography monitors
- Designed with 40% less dead space and superior oxygen delivery compared to blow-by designs
- Ergonomic design provides superior comfort and fit on both adult and pediatric faces
- Proprietary hydrophobic 0.2 micron filter protects monitor from biohazards
- Nafion® tubing for high humidity applications



Divided nasal version available

Comparable to standard orange and long term yellow reflective connectors



Nafion tubing for high humidity applications

END TIDAL SAMPLING LINE

ITEM	DESCRIPTION	LENGTH	PK	\$
4MSF5-13-25	ETCO ₂ Sampling Line with 15 mm Adapter	13'	25	\$239.95
4MSF5-6-25	ETCO ₂ Sampling Line with 15 mm Adapter	6'	25	\$216.95
4MSF8-13-25	ETCO ₂ Sampling Line	13'	25	\$152.95
4MSF8-6-25	ETCO ₂ Sampling Line	6'	25	\$152.95



DIVIDED ETCO₂ / O₂ NASAL CANNULA

- True ETCO₂ sampling and provided O₂ delivery without limitations at higher O₂ rates
- Soft tubing for patient comfort
- Divided design separates gases
- Salter Eyes maintains sampling even during a cannula occlusion



ITEM	DESCRIPTION	LENGTH	PK	\$
4950FT-7-7-25	Divided Cannula with Male Luer Connector	7'	25	\$49.95
4950FT-TG7-7-25	Divided Cannula with Threaded Grip Connector	7'	25	\$50.95

STAT ORANGE CANNULA

ITEM	DESCRIPTION	SIZE	TYPE	O ₂ X CO ₂	PK	\$
4MSF1-0-6-25	O ₂ /ETCO ₂ STAT Orange	Adult	Oral/Nasal	2' x 6'	25	\$181.95
4MSF1-7-13-25	O ₂ /ETCO ₂ STAT Orange	Adult	Oral/Nasal	7' x 13'	25	\$207.95
4MSF1-7-6-100	O ₂ /ETCO ₂ STAT Orange	Adult	Oral/Nasal	7' x 6'	100	\$756.95
4MSF1-L-7-6-25	O ₂ /ETCO ₂ STAT Orange	Adult	Oral/Nasal	7' x 6'	25	\$296.95
4MSF2-PED-7-6-25	O ₂ /ETCO ₂ STAT Orange	Pediatric	Oral/Nasal	7' x 6'	25	\$178.95
4MSF3-6-25	ETCO ₂ STAT Orange	Adult	Nasal	6'	25	\$167.95
4MSF3-7-6-25	O ₂ /ETCO ₂ STAT Orange	Adult	Nasal	7' x 6'	25	\$175.95



As Low As
\$1.99ea



VENTFLO ETCO₂ / O₂ SAMPLING CANNULA + REFLECTIVE CONNECTOR

- Dual nasal cannula simultaneously provides oxygen to both nostrils while obtaining CO₂ sampling during spontaneous breathing
- Designed to prevent mixing of fresh oxygen with CO₂
- Hydroscopic filter helps minimize moisture from being pulled into the monitor at the sampling port
- Innovative scoop design suitable for mouth-breathing patients
- 3-Channel tubing prevents kinking




~ Filter Integrated with Reflective Connector Helps Prevent Moisture Collection

~ Compatible with Microstream™, Capnostream™ and Oridion™ Capnography Monitors

~ Adult and Pediatric Sizes in 3 Tubing Lengths

~ DEHP Free



ETCO₂ / O₂ CANNULA + REFLECTIVE CONNECTOR

ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
5107F-SE	ETCO ₂ / O ₂ Nasal Cannula	Adult	7'	25	\$172.95
5110F-SE	ETCO ₂ / O ₂ Nasal Cannula	Adult	10'	25	\$190.95
5114F-SE	ETCO ₂ / O ₂ Nasal Cannula	Adult	14'	25	\$208.95
5207F-SE	ETCO ₂ / O ₂ Nasal Cannula	Pediatric	7'	25	\$218.95
5210F-SE	ETCO ₂ / O ₂ Nasal Cannula	Pediatric	10'	25	\$240.95
5214F-SE	ETCO ₂ / O ₂ Nasal Cannula	Pediatric	14'	25	\$265.95

ETCO₂ / O₂ CANNULA + REFLECTIVE CONNECTOR + SCOOP

ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
5707F-SE	ETCO ₂ / O ₂ Oral/Nasal Cannula	Adult	7'	25	\$203.95
5710F-SE	ETCO ₂ / O ₂ Oral/Nasal Cannula	Adult	10'	25	\$224.95
5714F-SE	ETCO ₂ / O ₂ Oral/Nasal Cannula	Adult	14'	25	\$246.95

As Low As
\$6.91ea



Strategically positioned holes deliver O₂ and collect CO₂ simultaneously through the nares

Soft-Ears™ material ensures superior patient comfort, even during prolonged use



As Low As
\$8.15ea



Fits-All O₂ connector is standard on all sizes and lengths


Scoop design effectively captures end tidal sampling from mouth breathing patients

RESPIRATORY

RESPIRATORY




Less than \$2.39 ea




ventlab
END TIDAL CO₂ + O₂ DELIVERY CANNULA

- Dual nasal cannula simultaneously provides oxygen to both nostrils while obtaining CO₂ sampling during spontaneous breathing
- Effectively captures undiluted end tidal gas sample
- Designed to prevent mixing of fresh oxygen with CO₂
- 3-Channel tubing prevents kinking
- Includes male/female adapter with luer lock connector
- Available with Fits-All Connector for O₂ or oxygen connections
- Available in three tubing lengths

ITEM	DESCRIPTION	TUBING LENGTH	PK	\$
4107	Adult ETCO ₂ Cannula	7'	25	\$59.95
4107F	Adult ETCO ₂ Cannula + Fits-All Connector	7'	25	\$62.95
4110	Adult ETCO ₂ Cannula	10'	25	\$69.95
4110F	Adult ETCO ₂ Cannula + Fits-All Connector	10'	25	\$73.95
4114	Adult ETCO ₂ Cannula	14'	25	\$79.95
4114F	Adult ETCO ₂ Cannula + Fits-All Connector	14'	25	\$82.95



Less than \$2.99 ea



ventlab
END TIDAL CO₂ + O₂ DELIVERY + SCOOP CANNULA

- Innovative design suitable for mouth breathing patients
- Design helps lessen attenuation of the end tidal sample by oxygen dilution
- Dual nasal cannula simultaneously provides oxygen to both nostrils while obtaining CO₂ sampling during spontaneous breathing
- Designed to prevent mixing of fresh oxygen with CO₂
- 3-Channel tubing prevents kinking
- Includes male/female adapter with luer lock connector
- Available in three tubing lengths

ITEM	DESCRIPTION	TUBING LENGTH	PK	\$
4707	Adult ETCO ₂ Cannula + Scoop	7'	25	\$74.95
4710	Adult ETCO ₂ Cannula + Scoop	10'	25	\$84.95
4714	Adult ETCO ₂ Cannula + Scoop	14'	25	\$94.95




SOFTECH® BI-FLO™ CANNULA

- Allows delivery of oxygen to both nasal prongs and sampling of expired gases from both nares simultaneously
- Available in three convenient adult configurations for use with any side stream CO₂
- Each cannula packaged with adapter
- Complete with gas sampling line

ITEM	OXYGEN TUBING	SAMPLING LINE	LUER CONNECTOR	PK	\$
8-3558-43	7'	3"	Female	25	\$96.95
8-3558-44	7'	7"	Male	25	\$96.95
8-3558-45	7'	7"	Large Female	25	\$96.95

Only \$3.87 ea




END TIDAL CO₂ + O₂ DELIVERY CANNULA

- Two end tidal CO₂ ports provides excellent sampling accuracy and superior waveforms
- Reliable sampling regardless of clogged sinus in one nostril, deviated septum or polyps
- Oxygen delivered at a lower rate, thus minimizing mucosal wall damage
- Secure luer lock connection for both male and female connectors
- Soft tip cannula provides superior patient comfort and nare positioning
- Non-kinking tubing ensures consistent and accurate delivery
- Universal oxygen connector standard on all versions

ITEM	OXYGEN TUBING	LUER LOCK	DESCRIPTION	PK	\$
8-3534-02	7'	Male, 2" CO ₂ Line	Adult ETCO ₂ Cannula	25	\$101.95
8-3534-04	7'	Female, Short CO ₂ Line	Adult ETCO ₂ Cannula	25	\$101.95
8-3534-06	7'	Male	Adult ETCO ₂ Cannula	25	\$101.95
8-3534-08	7'	Female	Adult ETCO ₂ Cannula	25	\$101.95
8-3534-10	10'	Male	Adult ETCO ₂ Cannula	25	\$106.95
8-3534-12	10'	Female	Adult ETCO ₂ Cannula	25	\$106.95

As Low As \$4.07 ea



CAPNOGRAPHY CANNULA

- Designed with minimal dead space in sampling line; resulting in a clear, accurate capnogram
- Comfort Soft Plus® tubing delivers premium patient comfort, reducing pressure ulcers and promoting patient wellness
- Oral/nasal design for mouth-breathers

ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
0912	CO ₂ /O ₂ Nasal Cannula	Adult	10'	25	\$246.95
0913	CO ₂ /O ₂ Oral/Nasal Cannula	Adult	10'	25	\$260.95
0914	CO ₂ /O ₂ Nasal Cannula	Adult	14'	25	\$258.95
0915	CO ₂ /O ₂ Oral/Nasal Cannula	Adult	14'	25	\$274.95
0916	CO ₂ /O ₂ Nasal Cannula	Pediatric	10'	25	\$235.95
0917	CO ₂ /O ₂ Oral/Nasal Cannula	Adult	7'	25	\$263.95
0918	CO ₂ /O ₂ Oral/Nasal Cannula with 2" Pigtail Connector	Adult	7'	25	\$260.95

As Low As \$9.44 ea



RESPIRATORY



DISCOUNT PRICING AVAILABLE!
CALL 800-237-5481

NOW ONLY \$44
FH603006

ventlab HUMIDIFIERS (HMEs)

- Bacterial/viral filters for use with PEEP valves and resuscitation devices
- Humidity moisture exchangers (HME) and filters (HMEF)
- Effectively provide heat and moisture delivery/recovery for ventilated and spontaneously breathing patients
- Compatible with leading breathing systems and come with or without filtration material
- 99.99% Filtration efficacy
- Accessory flex tubes and cuff available



HEAT MOISTURE EXCHANGER (HME)						
ITEM	FH603007 / FH603006 Port / No Port	FH603010	FH603012	FH603014	FH603021 / FH603023 Port / No Port	FH603028
DESCRIPTION	Round, Large	Swivel	Rectangular	Angled, Small	Straight, Small	Straight, Pediatric
VT RANGE (mL)	150 - 1500	250 - 1500	200 - 1500	150 - 1500	150 - 1500	150 - 1500
MOISTURE OUTPUT (mg H ₂ O/L) ±0.1g	33.4	34	32	33.2	33.7	34.0
WEIGHT (g)	17.8	37.1	36.4	16.8	14.7	13.6
ID / OD (mm)	22 ID x 22 OD / 15 ID	22 ID x 22 OD / 15 ID	15 ID x 22 OD / 15 ID	22 ID x 22 OD / 15 ID	22 ID x 22 OD / 15 ID	22 ID x 22 OD / 15 ID
DEAD SPACE (mL)	47	74	62	42.8	35.3	24.3
RESISTANCE (cm H ₂ O) 30 lpm	0.2	<1	<1	0.5	0.6	0.4
\$ PER CS OF 50	\$74.95 / \$74.95	\$99.95	\$74.95	\$74.95	\$74.95 / \$74.95	\$53.95

ventlab HEPA FILTER

- Effectively filters for COVID-19
- 22 mm OD/15 mm ID x 22 mm ID/15 mm OD
- Sampling port
- 99.9999% Filtration efficacy
- 62 mL Deadspace
- 250 mL – 1500 mL Tidal volume



NOW ONLY \$95
YOU SAVE \$75

HUDSON RCI BACTERIAL / VIRAL FILTER

- Bi-directional flow-low resistance
- Naturally hydrophobic filter media
- Specifications: bacterial filtration efficiency 99.999+%; viral-99.99+%
- Dead space 42ml
- Flow resistance: 1.26 cm H₂O at 60 l/min



ITEM	DESCRIPTION	PK	\$
FH603074	HEPA / COVID-19 Filter	50	\$169.95
3-5016-05	22mm ID x 15mm ID/22mm OD	50	\$106.95

ventlab FILTRATION & HUMIDIFICATION (HMEFs)

- Bacterial/viral filters for use with PEEP valves and resuscitation devices
- Humidity moisture exchangers (HME) and filters (HMEF)
- Effectively provide heat and moisture delivery/recovery for ventilated and spontaneously breathing patients
- Compatible with leading breathing systems and come with or without filtration material
- 99.99% Filtration efficacy
- Accessory flex tubes and cuff available



NOW ONLY \$46
FH603008

HEAT MOISTURE EXCHANGER + FILTRATION MEDIA (HMEF)						
ITEM	FH603008 / FH603005 Port / No Port	FH603009	FH603011	FH603013	FH603020 / FH603022 Port / No Port	FH603027
DESCRIPTION	Round, Large	Swivel	Rectangular	Angled, Small	Straight, Small	Straight, Pediatric
VT RANGE (mL)	150 - 1500	250 - 1500	200 - 1500	150 - 1500	150 - 1500	150 - 1500
MOISTURE OUTPUT (mg H ₂ O/L) ±0.1g	34.1	36	32.5	33.4	33.6	34.0
WEIGHT (g)	18.5	38.3	37.2	17.2	15.0	13.6
ID / OD (mm)	22 ID x 22 OD / 15 ID	15 ID x 22 OD / 15 ID	22 ID x 22 OD / 15 ID	22 ID x 22 OD / 15 ID	22 ID x 22 OD / 15 ID	22 ID x 22 OD / 15 ID
DEAD SPACE (mL)	47	74	62	42.8	35.5	24.3
RESISTANCE (cm H ₂ O) 30 lpm	0.8	<2	<1	1.9	2	2.4
B/V EFFICIENCY	>99.99%	>99.99%	>99.99%	>99.99%	>99.99%	>99.99%
\$ PER CS OF 50	\$75.95 / \$99.95	\$99.95	\$75.95	\$76.95	\$75.95 / \$75.95	\$53.95

Teleflex GIBECK® HUMID-FLO HEAT & MOISTURE EXCHANGER (HME)

- Has two operating modes: HME and AEROSOL, simply rotate blue collar to switch
- Allows delivery of medicated aerosol treatments without breaking circuits
- Eliminates caregiver exposure to de-pressurizing circuits
- Reduces both cross contamination potential and circuit manipulation



ITEM	DESCRIPTION	PK	\$
1-7201-12	Gibeck HME	25	\$235.95
1-7201-14	NebTee (22 mm O.D. x 18 mm I.D. x 22 mm I.D.)	30	\$171.95
1-7201-18	MDI Adapter	30	\$106.95



RESPIRATORY



Westmed

EXPANDABLE ADULT CIRCUIT

- High quality double limb tubing with exceptional flexibility
- Wide variety of accessory combinations; *more available – just call!*
- ISO-standard connections
- Plain or ported wyes
- HMEs, HMEFs, bacterial/viral filters and HEPA filter options
- Standard and swivel gas sampling ports, lines and filters
- Anesthesia face mask options
- ECG leads



CIRCUITS WITH 3L BREATHING BAGS

ITEM	DESCRIPTION	LENGTH	BAG VOL	PK	\$
9239	Exp. Circuit, GSL, GSE, Par-Wye, 2 Filters, Lg Mask	72"	3L	20	\$211.95
9677	Exp. Circuit, GSE, GSL, Par-Wye, Mini HME, 2 Filters, Lg Mask	96"	3L	10	\$68.95
8005	Exp. Circuit, Maxi HME, Lg Mask	96"	3L	20	\$320.95
8661	Exp. Circuit, GSE, GSL, Par-Wye, 2 Filters	96"	3L	20	\$182.95
9676	Exp. Circuit, GSE, GSL, Par-Wye, Mini HMEF, Lg Mask, ECG, Yankauer, Suction Tube	96"	3L	10	\$94.95
9344	Exp. Circuit, GSE, Par-Wye w/GSP, Mini HMEF, Filter, GSL M/F	96"	3L	20	\$180.95
9665	Exp. Circuit, GSL, GSE, Par-Wye, HMEF, Lg Mask	96"	3L	20	\$202.95
9261	Exp. Circuit, GSE, Par-Wye, Mini HMEF, Lg Mask, GSL	96"	3L	20	\$205.95
9110	Exp. Circuit, GSL, GSE, Par-Wye, Mini HMEF, Lg Mask	96"	3L	20	\$246.95
8143	Exp. Circuit, GSL, Angled HMEF, Swi-Wye, Lg Mask	108"	3L	20	\$227.95
9050	Exp. Circuit, GSE, Par-Wye w/port, Mini HMEF, GSL, Lg Mask	108"	3L	10	\$92.95
9611	Exp. Circuit, GSE, GSL, Par-Wye w/GSP, Mini HMEF, Lg Mask, Yankauer, Suction Tube	108"	3L	10	\$104.95
9362	Exp. Circuit, GSE, GSL, Par-Wye, Filter, Suction Tube	120"	3L	20	\$199.95

CIRCUITS WITH 2L BREATHING BAGS

ITEM	DESCRIPTION	LENGTH	BAG VOL	PK	\$
9574	Exp. Circuit, GSL, GSE, Par-Wye, 2 Filters	60"	2L	20	\$144.95
9342	Exp. Circuit, GSE, Par-Wye w/GSP, Filter	96"	2L	20	\$136.95
9671	Exp. Circuit, GSE, GSL, Par-Wye, Maxi HMEF	96"	2L	20	\$179.95
9129	Exp. Circuit, GSL, Mini HMEF	96"	2L	20	\$177.95
9246	Exp. Circuit, GSL, Par-Wye w/GSP, Filter, Lg Mask	96"	2L	20	\$259.95
9623	Exp. Circuit, GSE, GSL, Par-Wye, 2 Filters, Lg Mask, Yankauer, Suction Tube	108"	2L	10	\$123.95
9063	Exp. Circuit, GSL, GSE, Par-Wye, 2 Filters	120"	2L	20	\$252.95
9575	Exp. Circuit, GSE, GSL, Par-Wye, 2 Filters	120"	2L	20	\$157.95
9203	Exp. Circuit, GSE, Par-Wye	144"	2L	20	\$159.95
9636	Exp. Circuit, GSE, GSL, Par-Wye, 2 Filters, Lg Mask, Yankauer, Suction Tube	200"	2L	10	\$155.95



ventlab

BREATHING CIRCUIT

- Use as a conduit for gases between a Respironics® positive pressure generator and a patient during non-invasive ventilation
- 72" smooth bore tubing with pressure monitoring line and dual ported adapter
- 24" extension tubing with T-assembly for dual pressure monitoring



ITEM	DESCRIPTION	PK	\$
RC1000	Breathing Circuit with 22 mm OD Adapter	10	\$221.95
RC1000F	Breathing Circuit with 22 mm OD Adapter and Filter	10	\$221.95

FREE SHIPPING
for ONLINE ORDERS
of \$500+
www.bay-medical.com

Some restrictions apply.



Westmed

MR CONDITIONAL TRANSPORT CIRCUIT


- Complete, lightweight, ready-to-use, non-metallic MR conditional transport circuit in a single small package
- Compatible with some portable transport ventilators that have a passive exhalation valve
- Features a minimal deadspace, non-rebreathing valve with a swivel patient connector, adjustable PEEP valve and a proximal monitoring port that can be used for airway pressure
- Circuit expands up to 8 feet to allow convenient ventilator positioning away from the patient or outside of the MRI magnet
- Single-piece patient valve with integral adjustable PEEP capability
- Monitoring port with tethered cap allows connection of external pressure gauge or sidestream capnograph for expanded monitoring capability



ITEM	DESCRIPTION	SIZE	LENGTH	PK	\$
9317	MR Conditional Circuit, Patient Valve, PEEP	Adult	96"	10	\$210.95
9318	MR Conditional Circuit, Patient Valve, PEEP, CO ₂ Easy	Adult	96"	10	\$259.95

CIRCUIT EXTENDERS

ITEM	DESCRIPTION	SIZE	PK	\$
7-5132	Flex Extender	5.2"	50	\$74.95
S6411P	Flex Extender	5.2" x 15 mm OD x 15 mm ID	50	\$68.95
S6419P	Flex Extender	5.2" x 15 mm OD x 22 mm x 15 mm ID	50	\$64.95
S6420P	Flex Extender with Capd GSP	5.2" x 15 mm OD x 22 mm x 15 mm ID	50	\$112.95



Westmed

UNI-LIM ADULT CO-AXIAL CIRCUIT

- The optimal choice in co-axial circuits
- Streamlined tube-within-a-tube design that eliminates twisting between the inspiratory and expiratory limbs
- The internal inspiratory limb is warmed by the outside expiratory limb
- Provides natural heat and moisture conservation properties comparable to that of a standard circle system
- Designed to function as a standard circle anesthesia breathing system
- Eliminates bulky connectors at the patient interface
- Engineered with a secure-fit inspiratory limb to prevent patient-end disconnects
- No "floating" internal tubes
- ISO-standard connections for secure leak-resistant attachments



ITEM	DESCRIPTION	LENGTH	BAG VOL	PK	\$
9432	Uni-Lim Circuit, GSE, Maxi HMEF	96"	None	20	\$238.95
9061	Uni-Lim Circuit, GSE, GSL, 2 Filters, Adult Mask	60"	2L	20	\$304.95
9247	Uni-Lim Circuit, Elbow, HEPA, GSL, Yankauer, Suction Tube, ECG	72"	2L	20	\$426.95
9128	Uni-Lim Circuit, GSE, GSL, Mini HMEF	72"	2L	20	\$195.95
9615	Uni-Lim Circuit, GSE, Mini, HMEF, Lg Mask, Suction Tube, Yankauer, GSL	60"	3L	10	\$104.95
9566	Uni-Lim Circuit, Elbow, Maxi HMEF, IGSL, GSL	120"	3L	20	\$317.95
9106	Uni-Lim Circuit, GSE, GSL, 2 Filters, Lg Mask	60"	3L	20	\$233.95
9634	Uni-Lim Circuit, Elbow, Maxi HMEF, IGSL, GSL, Lg Mask	60"	3L	20	\$292.95
9170	Uni-Lim Circuit, GSE, GSL, Mini HMEF, Lg Mask	60"	3L	20	\$247.95
9153	Uni-Lim Circuit, GSE, GSL, HMEF, Mask	60"	3L	20	\$255.95
9565	Uni-Lim Circuit, Elbow, HMEF, IGSL	60"	3L	20	\$268.95
9096	Uni-Lim Circuit, GSE, Filter	120"	3L	20	\$196.95

RESPIRATORY



HIGH CONCENTRATION NON-REBREATHER & PARTIAL REBREATHER MASKS

- Elongated face mask design
- Safety vent option available
- Elastic head strap and over-the-ear styles
- Adult & pediatric face masks



ITEM	DESCRIPTION	INCLUDES	SIZE	PK	\$
8130-9-25	Non-Rebreather Mask	9' Tubing, Safety Vent	Adult	25	\$32.95
1130TG-7-50	Non-Rebreather Mask	7' Tubing, Threaded Grip, Safety Vent	Pediatric	50	\$126.95
8030-7-50	Non-Rebreather Mask	7' Tubing, Safety Vent, Cannula Style Strap	Adult	50	\$48.95
8040-7-50	Non-Rebreather Mask	7' Tubing, Cannula Style Securement	Adult	50	\$64.95
8130TG-7-50	Non-Rebreather Mask	7' Tubing, Threaded Grip, Safety Vent, Elastic Strap	Adult	50	\$64.95
1102-7-50	Partial Rebreather Mask	7' Tubing	Adult	50	\$35.95
1127-7-50	Partial Rebreather Mask	7' Tubing	Pediatric	50	\$99.95



SOFT-TO-THE-TOUCH OXYGEN CANNULA

- Face piece, nares and entire length of tubing is ultra soft, delivering exceptional patient comfort
- Soft, pliable, stretchable tubing
- Slide bolo prevents twisting and turning



ITEM	DESCRIPTION	LENGTH	PK	\$
1650FT-4-50	Adult O ₂ Cannula	4'	50	\$22.95
1650FT-7-50	Adult O ₂ Cannula	7'	50	\$25.95
1650FT-14-50	Adult O ₂ Cannula	14'	50	\$37.95
1650FT-25-50	Adult O ₂ Cannula	25'	25	\$39.95
1650FTTG-7-50	Adult O ₂ Cannula & Threaded Grip	7'	50	\$22.95
1650FTTG-14-50	Adult O ₂ Cannula & Threaded Grip	14'	50	\$25.95
1650FTTG-25-50	Adult O ₂ Cannula & Threaded Grip	25'	25	\$37.95



REUSABLE MASK HARNESS

- Silicone
- Autoclavable up to 134°C!



ITEM	DESCRIPTION	PK	\$
4-1050-12	Pediatric	10	\$70.95
4-1050-13	Adult	10	\$70.95



MASK HARNESS

- Reinforced holes on adult size for durability
- Center vents for comfort and fit at back of head
- Flexible neoprene
- Molded, not stamped, individually packaged



ITEM	DESCRIPTION	PK	\$
4-1054-23	Adult Mask Harness	1	\$13.95
4-1054-22	Child Mask Harness	1	\$13.95
4-1054-21	Infant Mask Harness	1	\$13.95
4-1054-23-10	Adult Mask Harness	10	\$98.95
4-1054-22-10	Child Mask Harness	10	\$98.95
4-1054-21-10	Infant Mask Harness	10	\$98.95




MASK HARNESS

- Light blue molded head harness
- Equal distribution of pressure over large area
- Designed to ensure proper fit



ITEM	DESCRIPTION	PK	\$
4-1054-01	Adult Blue Mask Harness	1	\$5.95



FACE MASK

- Inflatable cuff ensures seal and patient comfort
- Removable color-coded hook ring for hand-held procedures
- Adjustable cushion assures secure fit and comfort
- Transparent mask enables patient monitoring



ITEM	DESCRIPTION	HOOK RING COLOR	FITTING	PK	\$
4-0110-00	Neonate	Gray	15mm Male	30	\$74.95
4-0110-01	Infant	Lt Blue	15mm Male	30	\$74.95
4-0110-02	Child	Yellow	22mm Female	30	\$74.95
4-0110-03	Small Adult	Green	22mm Female	30	\$74.95
4-0110-04	Medium Adult	Red	22mm Female	30	\$74.95
4-0110-05	Large Adult	Dk Blue	22mm Female	30	\$74.95



RESPIRATORY

RESPIRATORY



NebuTech's Effective Drug Delivery & Patient Adaptability Can:

- Lower number of patient treatments required
- Reduce overall patient stay
- Help lower patient re-admission rates
- Promote patient compliance with in-home treatment



NEBUTECH® NEBULIZER

- Delivers front-loaded bolus dosing during critical first 1/3 of each breath
- Provides drug deposition deep into the lungs
- High quality aerosol – 80% of particles in the respirable range (<5 microns)
- Improved patient outcomes translate into financial benefits
- Quick yet effective treatment – 50% shorter treatment times than AeroEclipse II
- 82% of treatments completed in less than 3 minutes
- 45° angle for maximizing patient comfort and treatment compliance
- Patient-friendly accessories provide flexible and improved respiratory treatment
- Filter provides RT protection from exhaled medications
- Spans the full continuum of patient care: pediatrics through elderly



ITEM	DESCRIPTION	PK	\$
8992TG-7-50	NebuTech, Mask, 7' Tubing, TG Connector, B/V Adapter, Expiratory Filter	50	\$268.95
8990-7-50	NebuTech, Mouthpiece, 7' Tubing, B/V Adapter, Expiratory Filter	50	\$329.95
8990TG-7-50	NebuTech, Mouthpiece, 7' Tubing, TG Connector, B/V Adapter, Expiratory Filter	50	\$329.95
8992-7-50	NebuTech, Mask, 7' Tubing, B/V Adapter, Expiratory Filter	50	\$268.95



T-PIECE NEBULIZER

- Reliable and cost effective aerosol therapy
- Small volume jet nebulization
- Designed with anti-drool "T" for extended, uncontaminated use
- Ready to use with T-piece, mouthpiece, reservoir tube and supply tubing
- Easy to use; supports patient compliance
- Nebulizes faster than other brands of T-piece devices
- Designed for patient comfort during treatment



ITEM	DESCRIPTION	PK	\$
8900-7-50	Nebulizer with T, Mouthpiece, 6" Reservoir Tube, 7' Tubing	50	\$38.95
8900TG-7-50	Nebulizer with T, Mouthpiece, 6" Reservoir Tube, 7' Tubing, Thread Grip	50	\$38.95



HOPE™ NEBULIZER

- Refillable, cost effective, high output, continuous medication nebulizer
- Flow design and 200ml reservoir allow for up to 8 hours of nebulized medication delivery at 25ml/hr



ITEM	DESCRIPTION	PK	\$
3-1132-00	High Output Nebulizer	25	\$328.95



BREATHE RITE™ VALVED HOLDING CHAMBER

- Anti-static chamber enables effective medication delivery
- Designed to promote greater drug volume output
- Inspiratory aperture ensures proper inspiratory rate, important for complete delivery deposition
- Chamber design eliminates the need for a flow signal since whistling signals are shown to confuse users
- Collapsible version makes it very portable
- Face mask (all detachable) options available



ITEM	DESCRIPTION	MASK SIZE	PK	\$
BR200	Collapsible VHC	-	20	\$105.95
BR300	Rigid VHC	-	20	\$105.95
BR302-I	Rigid VHC with Mask	Infant	20	\$160.95
BR303-I	Rigid VHC with Mask	Small Child	20	\$160.95
BR304-I	Rigid VHC with Mask	Child	20	\$160.95
BR305-I	Rigid VHC with Mask	Adult	20	\$160.95



BUBBLE HUMIDIFIERS

- Dependable, cost effective humidification
- Compatible with all oxygen sources: concentrators, liquid tanks, wall outlets
- 360° diffusion port offers quiet operation with improved hydration and less back-pressure
- 3 and 6 PSI models
- High flow unit accommodates 6 to 15 LPM with audible pressure relief



ITEM	DESCRIPTION	PK	\$
7100-0-50	350 cc Bubble Humidifier, 3 PSI	50	\$80.95
7600-0-50	350 cc Bubble Humidifier, 6 PSI	50	\$80.95
7900-0-25	350 cc High Flow Bubble Humidifier	50	\$87.95



HIGH FLOW O₂ DELIVERY WITH HUMIDIFICATION

- 6-15 LMP
- Humidifier features 6 PSI (414 mbar) pressure relief valve
- Curved nares made of ultra-soft material for patient comfort
- 3-Channel tubing helps prevent kinking
- Green supply tubing standard
- Cannula features Fits-All Connector



ITEM	DESCRIPTION	PK	\$
6900	Bubble Humidifier	10	\$27.95
6907	7' High Flow O ₂ Cannula & Bubble Humidifier	10	\$47.95
6915	15' High Flow O ₂ Cannula & Bubble Humidifier	10	\$61.95



RESPIRATORY

RESPIRATORY



JET VENTILATOR

- Can be used in conjunction with needle cricothyroidotomy to ventilate patients
- Ideal for urgent care facilities, operating rooms, crash carts and EMTs
- Includes: Jet Ventilator, regulator, 6' high pressure hose and 4' small bore tubing

ITEM	DESCRIPTION	PK	\$
9-0004-00	Manual Jet Ventilator Kit	1	\$405.95

ACCESSORIES

ITEM	DESCRIPTION	PK	\$
9-0004-14	Small Bore Tubing, 4'	1	\$23.95
9-0004-16	High Pressure Hose, 6'	1	\$42.95
9-0004-05	Replacement Regulator	1	\$200.95




TRANSTRACHEAL CATHETER

- For use with jet ventilators
- 15mm and Luer lock connection
- Foam neck tape with Velcro lock
- Available in 13G and 14G catheters



ITEM	DESCRIPTION	SIZE	PK	\$
1-1804-02	Transtacheal Catheter	14G Child	1	\$111.95
1-1804-03	Transtacheal Catheter	13G Adult	1	\$111.95



JET VENT CATHETERS

- Standard and single lumen catheters available
- Allows a safe and atraumatic transtacheal access for oxygenation/ventilation with a manual jet or an automatic high frequency jet ventilator
- Three gauges available
- Does not include syringe



ITEM	DESCRIPTION	SIZE	PK	\$
9-0006-13	Jet Vent Catheter	13G-Adult	5	\$268.95
9-0006-14	Jet Vent Catheter	14G-Child	5	\$268.95
9-0006-16	Jet Vent Catheter	16G-Infant	5	\$268.95

9-0003-55



JET VENTILATOR


- A ventilating system for rigid-wall bronchoscopy laryngeal micro-surgery and emergency use
- Includes: 5' high-pressure tubing with 2 DISS fittings, flow control valve, thumb control valve, pressure gauge, and 20" clear plastic tubing with Luer lock fitting for attachment of adapters and tubes
- Custom tubing sizes available; please call for information

ITEM	DESCRIPTION	PK	\$
9-0003-25	Jet Ventilator System	1	\$460.95
9-0003-55	MRI-Safe Jet Ventilator System	1	\$1,462.95



EYEGARD™

- Clear plastic covers protect eyes during surgery and other medical procedures
- Safely adheres directly to patient skin
- Tab design means no more hassles with sticky tape and gloved hands
- Packaged in pairs



ITEM	DESCRIPTION	PK	\$
9-0210-50	EyeGard Eye Protector	50 Pr	\$47.95

9-0003-25



ACCESSORIES

ITEM	DESCRIPTION	SIZE	PK	\$
9-0003-27	Copper Tube, Infant	10cm	1	\$35.95
9-0003-28	Copper Tube, Child	15cm	1	\$35.95
9-0003-29	Copper Tube, Adult	20cm	1	\$35.95
9-0003-30	Copper Tube, Large Adult	30cm	1	\$48.95
9-0003-26	Jet Ventilator Bronchoscope Adapter	-	1	\$103.95
9-0003-31	Jet Ventilator Adapter 15mm F x Luer F	-	1	\$63.95
8-1400-08	High Pressure Hose	8'	1	\$94.95



RESPIRATORY

RESPIRATORY



7-1201-50



7-1201-80



7-1201-60



7-1201-70



7-1201-10



7-1201-40



7-1201-20



7-1201-30



CO₂ ABSORBENT

- Designed to eliminate the potential for production of compound A, carbon monoxide and excessive heat
- Features a permanent color change (off-white to violet) that will not revert over time; helps prevent inadvertent use of exhausted product
- Does not contain NaOH or KOH
- Suitable for low flow applications

ITEM	DESCRIPTION	PK	\$
7-1201-50	1 kg - Standard Systems with Cylindrical Canisters	12	\$185.95
7-1201-60	1.6 L Litholyme Refill	12	\$185.95
7-1201-70	GE® Multi-Style Cartridge - for Avance®, Aespire®, Aisys® Systems	6	\$165.95
7-1201-80	Dräger® Style Cartridge - for Apollo®, Pallas®, Primus® Systems	6	\$140.95



SODALIME CO₂ ABSORBENT

- Fully compliant U.S. Pharmacopoeia / National Formulary sodalime
- For anesthesiology, metabolators, oxygen therapy
- Prepack cartridge fits easily into Ohmeda and Dräger anesthesia machines
- Canister pack is in airtight plastic and foil lined bag



ITEM	DESCRIPTION	PK	\$
7-1200-40	Pre-pack Cartridge - 2.5 lbs.	12	\$256.95*
7-1200-30	Canister-pack - 3 lbs.	12	\$268.95*



CO₂ ABSORBENT


- Contains ethyl violet, which acts as a color indicator when absorbent is nearing exhaustion
- No KOH - minimum anesthetic agent degradation compared to other brands
- Low dust, high surface area and graded particle size
- Low risk of carbon monoxide formation due to good resistance of dry gas desiccation
- Low bulk density - less weight and less waste

ITEM	DESCRIPTION	PK	\$
7-1201-10	1 kg - Standard Systems with Cylindrical Canisters	12	\$117.95
7-1201-20	1.6 L Carbolime Refill	12	\$79.95
7-1201-30	GE® Multi-Style Cartridge - for Avance®, Aespire®, Aisys® Systems	6	\$106.95
7-1201-40	Dräger® Style Cartridge - for Apollo®, Pallas®, Primus® Systems	6	\$160.95



INFU-SURG® PRESSURE INFUSION BAG

- Pressure relief valve prevents over inflation (330 mmHg pressure relief)
- Large, oval shaped bulb allows for quick and easy inflation of bladder
- Single-handed inflation and deflation design makes it easy to use and requires minimal training
- Suitable for use with external inflation sources
- Color-coded gauge makes for accurate pressure monitoring (0-300 mmHg)
- Three-way stopcock ensures precise control of pressure
- Loads quickly and easily



Hook Style for 4030



Thread woven loop through infusion bag and hang on IV pole



ITEM	DESCRIPTION	SIZE	PK	\$
4005H	Infu-Surg Pressure Infuser Bag	500cc	5	\$65.95
4010H	Infu-Surg Pressure Infuser Bag	1000cc	5	\$71.95
4030	Infu-Surg Pressure Infuser Bag	3000cc	5	\$120.95



RESPIRATORY

FREE SHIPPING
for ONLINE ORDERS
of \$500+
www.bay-medical.com

Some restrictions apply.




MAGILL FORCEPS

- Reusable
- Wide variety of styles for practically every application
- Manufactured with surgical stainless steel

ITEM	DESCRIPTION	SIZE	PK	\$
9-3016-02	Magill Forceps, Open Tip, Child	4.5"	1	\$17.95
9-3016-03	Magill Forceps, Open Tip, Adult	5.5"	1	\$17.95
9-3018-01	Magill Forceps, Closed Tip, Infant	3"	1	\$17.95
9-3018-02	Magill Forceps, Closed Tip, Child	4"	1	\$17.95
9-3018-03	Magill Forceps, Closed Tip, Adult	5"	1	\$17.95




Closed Tip



Open Tip

INTUBATION




smiths medical ORIGINAL ESCHMANN ENDOTRACHEAL TUBE INTRODUCER (BOUGIE)

- Designed for changing or introducing full range of ET tubes
- Braided polyester for proper stiffness
- Coudé angled tip facilitates insertion in adults
- Individually packaged non sterile



Item	DESCRIPTION	Fr	size	PK	\$
9-0208-12	Pediatric, Straight Tip	10Fr	70cm	1	\$111.95
9-0208-13	Adult, Straight Tip	15Fr	70cm	1	\$111.95
9-0208-60	Adult, Angled Tip	15Fr	60cm	1	\$106.95



SunMed MALLEABLE INTRODUCER (BOUGIE)

- Combines benefits of an introducer and stylette into one unit
- Internal wire provides user the ability to pre-bend the stylette
- Graduation marks every 5cm
- Can be used with ET tubes, rigid bronchoscope, laryngeal type masks and double lumen endobronchial tubes



Item	DESCRIPTION	Fr	size	PK	\$
9-0213-82	Pediatric Malleable	10Fr	70cm	10	\$72.95
9-0213-92	Adult Malleable	15Fr	70cm	10	\$72.95

SunMed PORTED ENDOTRACHEAL INTRODUCER

- Provide oxygenation during intubation or remove light secretions
- Suction and extract tracheal gas samples for analysis
- Calibrated to ensure accurate distance insertion
- Fits in 6mm to 11mm tubes




Item	Description	Size	PK	\$
9-0212-80	Adult, Coudé Tip Ported Bougie	15Fr x 60cm	10	\$76.95

Now Only \$7.70ea

SunMed

ORIGINAL "BLUE" INTRODUCER (BOUGIE)

- Provides proper stiffness for ease of insertion
- Calibrated for easy observation of insertion
- Coudé tip facilitates insertion in adults
- Pediatric size has calibration marks on both ends allowing use as straight or Coudé tip introducer
- Individually wrapped



ITEM	DESCRIPTION	Fr	SIZE	Pk	\$
9-0211-70	Pediatric, Coudé/Straight Tip	10Fr	70cm	10	\$70.95
9-0212-72	Adult, Straight Tip	15Fr	70cm	10	\$70.95
9-0212-70	Adult Coudé Tip	15Fr	70cm	10	\$70.95

SunMed

BOUGIE-TO-GO™ ENDOTRACHEAL INTRODUCER

- Unwraps without distorting introducer shape; keeps tip in correct insertion position
- Low density polyethylene provides proper stiffness for ease of insertion
- Calibrated to ensure accurate distance insertion
- Fits in 6mm to 11mm tubes



Item	Description	Size	PK	\$
9-0212-82	Adult, Coudé Tip Bougie-To-Go	15Fr x 60cm	10	\$70.95

Conveniently packaged to fit in paramedic bags, pockets and compact spaces


Now Only \$7.10ea



INTUBATION

INTUBATION

INTUBATION



SunMed SUNGLIDE™ STYLETTE

- Malleable stylette - may be adjusted to any shape
- Stop enables the stylette tip to not protrude beyond the end of the ET tube
- Reduced surface tension removes need for lubricants
- Color coded stoppers
- Easy open package



Item	DESCRIPTION	Fr	FITS ET TUBE	LENGTH	PK	\$
9-0209-71	Neonatal, Blue Stopper	6Fr	2.0-3.5mm	13"	25	\$87.95
9-0209-72	Child, Gray Stopper	8Fr	4.0-7.0mm	15.5"	25	\$87.95
9-0209-73	Adult, Green Stopper	10Fr	7.0-10.0mm	18"	25	\$87.95

FREE SHIPPING
for ONLINE ORDERS
of **\$500+**
www.bay-medical.com

Some restrictions apply.



Instrumentation Industries, Inc.

ENDOTRACHEAL TUBE CHANGERS

- Simple, rapid technique for changing endotracheal tubes in less than two minutes
- Graduated marks for easy, accurate placement of new tube
- Available for all sizes of ET Tubes
- Size is clearly marked and easy to read
- Polyethylene



ITEM	DESCRIPTION	IMPRINT	FITS ET TUBE	LENGTH	PK	\$
9-0202-25	Endo Tube Changer	Black	2.5mm	15.25"	10	\$118.95
9-0202-30	Endo Tube Changer	Blue	3.0mm	16.5"	10	\$118.95
9-0202-40	Endo Tube Changer	Magenta	4.0mm	19.5"	10	\$118.95
9-0202-50	Endo Tube Changer	Blue	5.0mm	23"	10	\$118.95
9-0202-55	Endo Tube Changer	Black	5.5mm	25"	10	\$118.95
9-0202-60	Endo Tube Changer	Magenta	6.0mm	26"	10	\$118.95
9-0202-65	Endo Tube Changer	Green	6.5mm	26.5"	10	\$118.95
9-0202-70	Endo Tube Changer	Blue	7.0mm	27.5"	10	\$118.95
9-0202-99	Endo Tube Changer	Black	7.5-10.0mm	29.25"	10	\$118.95
9-0202-00	Complete set, one of each size				9	\$118.95



SunMed SUNSLIP SHEATHED STYLETTE

- Sealed in a PVC sheath
- Malleable aluminum



ITEM	DESCRIPTION	Fr	SIZE	PK	\$
9-0203-06	Infant SunSlip Stylette	6Fr	25cm	20	\$60.95
9-0203-10	Child SunSlip Stylette	10Fr	28cm	20	\$60.95
9-0203-14	Adult SunSlip Stylette	14Fr	32cm	20	\$60.95



RUSCH SLICK® STYLETTE ET TUBE GUIDE

- Facilitates first attempt intubation with minimal risk and trauma
- Easy open package



ITEM	DESCRIPTION	ET TUBE	COLOR	PK	\$
9-0209-51	Pediatric/Neonatal	2.0-3.5mm	Blue	25	\$93.95
9-0209-52	Small Adult/Child	4.0-7.0mm	Gray	25	\$94.95
9-0209-53	Adult	7.0-10mm	Green	25	\$94.95



SunMed SOLID COPPER STYLET

- Malleable
- Adjustable stop



ITEM	DESCRIPTION	LENGTH	PK	\$
9-0203-00	Solid Copper Stylet	14 3/4"	1	\$12.95


SunMed HALLOW COPPER STYLETS

- Hollow copper tube with suction port



ITEM	DESCRIPTION	LENGTH	PK	\$
9-0206-01	Hollow Copper Stylet with Suction	16 1/2"	1	\$27.95

INTUBATION



NOW ONLY \$60 FOR EACH

SunMed
MALLEABLE FIBER OPTIC INTUBATING STYLETTE

- Battery operated - 2 "AA"
- Supplied with adjustable stop
- Reusable

ITEM	DESCRIPTION	LENGTH	PK	\$
9-0204-18	Malleable F/O Stylette, Adult	18"	1	\$84.95
9-0204-14	Malleable F/O Stylette, Pediatric	14"	1	\$84.95
9-0210-05	Replacement Stops	-	5	\$6.95
5-0240-59	Xenon Fiber Optic 3.5mm Lamp	-	1	\$19.95



SunMed
LIGHT GUIDED STYLETTE

- Aids in proper placement of ET tubes
- Effectively transilluminates soft tissue of the neck
- Provides minimal neck flexion in trauma cases
- Aids in blind intubation
- 4 year shelf life



ITEM	DESCRIPTION	LENGTH	PK	\$
9-0210-01	Light Guided Stylette, Adult	15.75"	1	\$36.95



AINCA
VIDEO STYLETTE

- Inexpensive, lightweight, and reusable
- Allows practitioner to obtain a direct view from the tip of the ET tube as it progresses during intubation
- 5.8mm diameter and malleable along its entire length
- Powerful camera assembly with ultra-bright LED illumination at the tip delivers a crisp 2.4" color image on lightweight rechargeable LCD monitor
- Waterproof and can be cleaned by scrub/wash/wipe, high-level disinfected in 2% Glutaraldehyde, or sterilized in Sterrad® gas plasma hydrogen peroxide

ITEM	DESCRIPTION	LENGTH	FITS TUBE	PK	\$
9-0208-90	Video Stylette	12.3" / 32.0 cm	7.0 - 7.5mm	1	\$420.95

Powerful camera assemble delivers crisp 2.4" color image on lightweight rechargeable monitor



Ultra-bright LED illumination at the tip delivers crisp color image on LCD monitor

SURGICAL



MULTI-COMPARTMENT TRAYS

- Multi-compartment trays can be used inside drawers or directly on work surface
- Sturdy trays can be used to transport instruments and surgical components
- rounded corners for easy cleaning and smooth, non-scratching surface
- Constructed of high impact white polystyrene

ITEM	DESCRIPTION	SIZE	PK	\$
C8-2020-03	4 Compartment Tray	17" x 2" x 13"	1	\$56.95
C8-2020-04	DISC	14" x 3" x 10"	1	\$24.95
C8-2020-06	DISC	12" x 2" x 16"	1	\$28.95
C8-2020-01	7 Compartment Tray	17" x 2" x 13"	1	\$56.95
C8-2020-02	12 Compartment Tray	17" x 2" x 13"	1	\$56.95



CUSTOMIZABLE BINS

- Adjustable dividers allow for up to 16 bins
- Create bins of varying sizes for storage of small items and loose supplies
- Vertical orientation and small footprint means more space for working
- Constructed of durable white PVC

ITEM	DESCRIPTION	SIZE	PK	\$
C8-1010-02	2 Shelves, 14 Adjustable Bins	24" x 20" x 10"	1	\$249.95
C8-1010-03	16 Adjustable Bins	12" x 12" x 7"	1	\$169.95
C8-1010-01	1 Shelf, 10 Bins	24" x 16" x 11"	1	\$216.95





PRE-CUT LABELS, LABEL TAPE

- Meet ASTM Standards with proper drug labeling
Both labels and tape are color-coded by drug category
Show expiration date, expiration time, and initials of healthcare worker
Packaged on rolls for easy use and dispensing
Dispensers available on page xxx
Customized labels available for a minimum 6 roll purchase; call for more information

Table with 4 columns: ITEM, DESCRIPTION, PK, \$
SEE BELOW Pre-cut Labels 1 Roll \$10.95
SEE BELOW Label Tape 1 Roll \$10.95

Shown actual size



Pre-cut Labels 1/2" x 1-1/2" 500" Roll



Label Tape 1/2" x 1-1/2" 500" Roll

Narcotics

Blue with black imprint

Table with 3 columns: Pre-cut Labels, Label Tape, Description. Lists various narcotic drugs and their dosages.

Narcotic Antagonists

Blue and white diagonal stripe with black imprint

Table with 3 columns: Pre-cut Labels, Label Tape, Description. Lists narcotic antagonist drugs.

Muscle Relaxants

Red fluorescent with black imprint

Table with 3 columns: Pre-cut Labels, Label Tape, Description. Lists muscle relaxant drugs.

Relaxant Antagonists

Red and white diagonal stripe with black imprint

Table with 3 columns: Pre-cut Labels, Label Tape, Description. Lists relaxant antagonist drugs.

Miscellaneous

White with black imprint

Table with 3 columns: Pre-cut Labels, Label Tape, Description. Lists miscellaneous drugs.

Anticholinergic Agents

Green with black imprint

Table with 3 columns: Pre-cut Labels, Label Tape, Description. Lists anticholinergic agents.

Vasopressors

Violet with black imprint

Table with 3 columns: Pre-cut Labels, Label Tape, Description. Lists vasopressor drugs.

Hypotensive Agents

Violet and white diagonal stripe with black imprint

Table with 3 columns: Pre-cut Labels, Label Tape, Description. Lists hypotensive agents.

Local Anesthetics

Gray with black imprint

Table with 3 columns: Pre-cut Labels, Label Tape, Description. Lists local anesthetic drugs.

Sedatives

Orange with black imprint

Table with 3 columns: Pre-cut Labels, Label Tape, Description. Lists sedative drugs.

Major Sedatives

Salmon with black imprint

Table with 3 columns: Pre-cut Labels, Label Tape, Description. Lists major sedative drugs.

CENTRAL STERILE

CENTRAL STERILE



STERILIZATION INTEGRATORS

- Black indicator clearly denotes sterilization pass/fail
- Chemical indicator for monitoring steam sterilization processes - steam, time, temperature
- For use with steam sterilizers operating in the 250°F-275°F temperature range
- During sterilization, a black bar migrates along a paper wick through fail and PaSS areas indicating whether the sterilization conditions were met
- PaSS result assures that conditions in the pack were sufficient to kill even the most resistant forms of micro-organisms
- Meets requirements for Class 5 indicators according to aaMI/ISO 11140-1:2005

ITEM	DESCRIPTION	PK	\$
C6-1001-10	Vapor Line Sterilization Integrator	100	\$20.95
C6-1001-25	Vapor Line Sterilization Integrator	250	\$30.95



TEST PACKS FOR PRE-VACUUM STERILIZERS

- Detects leaks and low temperature pressure conditions, wet steam, super-heated steam and non-condensable gases
- Lead and heavy metal free; non toxic and safe to use
- Lot number marked on packs and test sheets
- Precise and clear color change indicators
- Test sheet can be archived as a performance record; ink will not revert



ITEM	DESCRIPTION	PK	\$
C6-3025-01	Once-A-Day Vertose Test Pack	30	\$148.95



SELF-CONTAINED BIOLOGICAL INDICATORS

- Available for steam or eO gas sterilizers
- Process indicator on each vial
- Distinctive color change of the indicator media for immediate and accurate reading
- Incubation time: 24hrs/steam vials; 48hrs/eO gas vials
- Color-coding caps and media: brown/steam; red/eO gas
- Universal indicator incubator designed for use with steam or eO gas vials
- Incubator holds up to 13 vials for incubation plus one crushing cavity
- Recommended temperatures: 57°C/steam; 37°C/eO gas; adjustable with toggle switch

ITEM	DESCRIPTION	PK	\$
C6-2000-00	Bi-OK Steam Vials	20	\$63.95
C6-2000-01	Bi-OK Steam Vials	100	\$286.95
C6-2001-01	Bi-OK EO Gas Vials	100	\$255.95
C6-2010-00	Universal Biological Indicator Incubator	1	\$618.95



STERILIZATION TUBES

- Steam and EO gas indicators both inside and outside
- Steam indicators: yellow changes to gray-black (red indicators do not change)
- EO gas indicators: red changes to green (yellow indicators do not change)
- Use a heat sealer to securely close



ITEM	DESCRIPTION	SIZE	PK	\$
C6-3010-02	Chex-All Sterilization Tube	2" x 100'	10	\$110.95
C6-3010-03	Chex-All Sterilization Tube	3" x 100'	10	\$124.95
C6-3010-04	Chex-All Sterilization Tube	4" x 100'	10	\$153.95
C6-3010-06	Chex-All Sterilization Tube	6" x 100'	10	\$189.95
C6-3010-09	Chex-All Sterilization Tube	9" x 100'	10	\$204.95
C6-3010-12	Chex-All Sterilization Tube	12" x 100'	10	\$274.95



NON-WOVEN STERILIZATION WRAP

- Designed for general hospital use and all types of packaging needs
- Delivers an efficient microbial barrier for pack security
- Tightly anchored fiber structure makes it virtually lint-free
- Does not delaminate; can use either side
- Excellent liquid repellency provides an effective barrier against contamination of contents

ITEM	DESCRIPTION	SIZE	PK	\$
C6-3020-12	Mini Steri-Wrap Sheets	12" x 12"	1000	\$125.95
C6-3020-15	Extra Small Steri-Wrap Sheets	15" x 15"	1000	\$165.95
C6-3020-24	Small Steri-Wrap Sheets	24" x 24"	500	\$198.95
C6-3020-36	Medium Steri-Wrap Sheets	36" x 36"	200	\$186.95
C6-3020-48	Large Steri-Wrap Sheets	48" x 48"	100	\$168.95



STEAM AUTOCLAVE TAPE

- Hold packs together securely without the need for pins, staples or fasteners
- Pressure sensitive tape imprinted with process indicator steam reactive ink
- When exposed to steam, the white/cream ink changes to brown/black; it is unaffected by dry heat or light
- Enables immediate identification between processed and unprocessed packs
- Adheres to a variety of sterilization packaging, including fabric and non-woven wraps
- Readily accepts ball point pen or pencil writing
- Meets requirements for Class 1 Process Indicators according to AAMI/ISO 11140-1:2005

ITEM	DESCRIPTION	WIDTH	PK	\$
C6-1010-50	Strate-Line Autoclave Tape, 60 yds	0.75"	500	\$206.95
C6-1010-60	Strate-Line Autoclave Tape, 60 yds	1.0"	200	\$206.95



CORNER PROTECTORS

- No more torn or ripped wrappings
- Prevents sterilization wrap around corners and blunt edges of trays from tearing or ripping
- Constructed of durable silicone
- Compatible with steam, gas, ozone, and low temperature plasma sterilization
- Reusable
- One size fits most surgical trays

ITEM	DESCRIPTION	PK	\$
C4-9000-01	Corner Protectors	80	\$244.95



INSTANT SEALING AUTOCLAVE POUCHES

- Steam and EO gas indicators both inside and outside
- Delivers considerable cost savings since a separate indicator does not need to be inserted into pouch
- Integrated adhesive strip eliminates the need for heat sealing
- Designed to be tamper-proof
- Steam indicators: yellow changes to gray-black (red indicators do not change)
- EO gas indicators: red changes to green (yellow indicators do not change)



ITEM	DESCRIPTION	SIZE	PK	\$
C6-3000-08	Chex-All II Instant Sealing Pouch	3" x 8"	500	\$69.95
C6-3000-10	Chex-All II Instant Sealing Pouch	5" x 10"	500	\$96.95
C6-3000-12	Chex-All II Instant Sealing Pouch	7" x 12"	500	\$136.95
C6-3000-14	Chex-All II Instant Sealing Pouch	8" x 16"	250	\$132.95
C6-3000-15	Chex-All II Instant Sealing Pouch	5" x 15"	500	\$197.95
C6-3000-18	Chex-All II Instant Sealing Pouch	12" x 18"	250	\$191.95

Customer I.D.

Priority Code



Promote safety and protection for every patient, every clinician. Order today.



Patient Eye Protection
Page 39



Face Protection
Page 39

Proper Labeling for Controlled Substances
Pages 64-65

Fentanyl _____ mcg/mL	
Prep Date _____	Exp. Date _____
Initial _____	Exp. Time _____



Sterilization Supplies
Pages 66-67



Filters
Pages 46-47



Single Patient Use Laryngoscopes
Pages 4-9